NSK

Multi Function Ultrasonic Scaler

Varios 970

Varios 970

iPiezo engine

OPERATION MANUAL

Please read this Operation Manual carefully before use,
and file for future reference.

Rx Only

OM-E0619E 001

MADE IN JAPAN

◆ **Classifications of equipment**

- Type of protection against electric shock:
 - Class I equipment
- Degree of protection against electric shock:
 - Type BF applied part: 
- Method of sterilization or disinfection recommended by the manufacture:
 - See 14. Sterilization
- Degree of protection against ingress of water as detailed in the current edition of IEC 60529:
 - Foot Control: IPX1 (Protected against vertically falling water drops)
- Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:
 - EQUIPMENT not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Mode of operation:
 - Continuous operation

Intended to Use

This product is designed only for dental clinic/dental office use. This device, Varios 970 / Varios 970 Lux, an electronic ultrasonic scaler, is intended for use with an appropriate tip for the following use: Scaling, Perio, Implant Maintenance, Endodontic, Retrograde Endo, Restorative (for Minimal Intervention/ Finishing/Trimming/ Polishing/Caries of Dentin), Prosthetics (Condensation / Loosening / Plugging).

Rx Only Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

1. ⚠ Cautions for handling and operation


- Read these cautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	A hazard that could result in bodily injury or damage to the device if the safety instructions are not followed.
 CAUTION	A hazard that could result in light or moderate bodily injury or damage to the device if the safety instructions are not followed.
 NOTICE	General information needed to operate the device safely.

 **WARNING**

- TO PREVENT ELECTRIC SHOCK Do not unplug the power cord with wet hands.
- TO PREVENT ELECTRIC SHOCK Make sure not to get water on the Control Unit.
- TO PREVENT ELECTRIC SHOCK Do not touch the handpiece backend electrical connections.
- TO PREVENT ELECTRIC SHOCK Use an electrical outlet that is grounded.
- If you feel any abnormality such as vibration, heat generation, abnormal noise, etc., prior or during the use of the unit, stop using it immediately.
- Do not turn the Power Switch without reason; it might blow out a fuse.
- This product is Medical Electrical equipment Electromagnetic compatible (EMC).As described in the accompanying documentation.
- Portable and mobile RF communications equipment can affect Electrical Medical equipment. Do not use RF equipment in close proximity to the product.
- When installing the product, provide space of approximately 10cm around the Control Unit for easy access to the inlet and the Power Cord.
- Staring directly into LED may result in eye damage. Do not turn it to the eyes of the patient or operator.

- Be sure to always use with sufficient coolant water. Insufficient coolant water could cause a burn of the affected area or the overheating of the instrument.
- Be aware of the tip end position. Constantly move the tip so that it does not stagnate at one point. If the tip is stagnated at a single point, there will not be enough cooling, and this could result in overheating of the tooth or root canal, or excessive cutting of the tooth or root canal.
- USE ONLY NSK genuine tips when using NSK Varios Ultrasonic Scaler (Varios 970 or Varios 970 LUX) problems such as damage, failure and accident of Handpieces resulting from use of Non-NSK Tips are not included in the warranty. The following are the possible failure that could happen when using the Non-NSK Tips;
 - Vibration failure caused by using non conforming screws.
 - Patients accidental ingestion of broken tips.
 - Damage of thread ridge of handpiece.
- You must use the tip within the power range described on the Tip-Power Guide. If you use it out of the power range, the tip might break or damage an operative site.
- When operating the product always consider the safety of the patient.
- Use by medical professional, such as doctor or dental hygienist, is intended.
- Check the vibration outside the patient's oral cavity before use. If any abnormalities are found, stop using immediately and contact your dealer.
- Do not drop, hit, or give excessive shock on the Control Unit/Handpiece. This could lead to electric shock or product failure. In case of giving shock to an optic handpiece, illumination parts that are equipped with the handpiece could be damaged and drop out of the handpiece.
- To prevent possible tooth plane damage and handpiece overheating, always use with sufficient water.
- Do not sterilize by ultraviolet light. Handpiece could discolor.
- Sterilize the Tip (except for all diamond coated tip), Handpiece, Tip Holder, Tip Cover S and Tip Wrench by autoclaving. Wipe the Control Unit, AC Power Cord, Foot Control, Tip Card, Tip Power Guide, and Handpiece Cord including the cover with alcohol soaked cotton.
- If chemical, solvent or antiseptic solution is deposited on this product, immediately wipe it away. Discoloration or deformation may occur if left.
- Do not disassemble or alter the handpiece/Control Unit.
- Keep away from patients with cardiac pacemakers.
- Keep away from explosive substances and flammable materials. Do not use for patients anesthetized under laughter gas. (Nitrous Oxide)
- Use the Fuse of specified rating. (120V: T630mA 250V, 230V: T315mA 250V)
- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information.
- Electromagnetic Compatibility: The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of this product as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of this product.
- This product should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this product should be observed to verify normal operation in the configuration in which it will be used.
- If any water drops remain on the handpiece after autoclaving, wipe them off., Staining may result if left.
- Do not use water tubes of other companies. Doing so may cause a malfunction.
- Do not spray using ozonated water, high oxidation potential water (strongly acidic water, superacidic water, electrolyzed water, etc.), or hypochlorous acid water. Doing so may lead to clogging, corrosion of metal parts, etc.
- Grounding reliability can only be achieved when the equipment is connected to an equipment receptacle marked "Hospital Only" or "Hospital Grade".
- Do not apply excessive power to the Tip. It may damage the teeth because of the ultrasonic vibration.




CAUTION

- During operation, high frequency oscillations in the handpiece and handpiece cord may affect computer and LAN Noise may be heard during operation near a radio receiver.
- Be sure to turn off the Power Switch after use. Remove the power plug and water inside of the Control Unit before storage.
- Users are responsible for operational control, maintenance and inspection.
- Clean/sterilize the product immediately after using it. Then store it. Leaving it non-sterile might lead to failure.
- When cleaning, use ethanol, or ethanol with isopropyl alcohol to wipe the surface of the handpiece. Use of chemicals other than those may cause the handpiece to discolor, crack, etc.
- When you have not used the product for long time and use it again, check the operation before use.
- This product does not consider patient's age (except infants), gender, weight or nationality.
- No special training is required for this device.
- Applied parts for patient and/or operator are/is tip and Handpiece.
- Surface temperature of tip shall be more than 50 degree without using a tap water or bottle. To avoid this event, be sure to use a tap water or bottle.
- When using a file as a tip for root canal treatment (for cleaning), use it after confirming its characteristics and areas where it can be applied.
- Operators and all others in the area must wear eye protection and a mask when operating this product.
- Before use, be sure to always check that the power level is set to the proper settings.
- During vibration, do not grasp the tip nor make contact with the tip for an extended time.
- Do not use tips that are damaged, bent, corroded, or show other imperfections. The tip could break or fail to vibrate.
- Do not apply more pressure than required to the tip end. This could cause the tip to break or warp. Do not twist or strain the tip, and do not use it as a replacement for surgical tools. This could cause the tip to deform or break.
- If the tip cracks or breaks during use, immediately stop operation, and replace with a new tip. Operation with a damaged tip could injure the gums, periodontal tissue, or other locations.
- The diamond coating on the tip surface may peel off depending on the tooth condition or root canal condition. If the diamond coating peels off, the cutting efficiency becomes less, resulting in reduced cutting performance. If you sense this type of condition, replace with a new tip.
- All sterilizable components of this device are delivered in a non-sterile condition and must be autoclaved prior to use. Refer to "5.Sterilization before first time use".
- The diamond coated tips are "Single use". Do not reuse.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

NOTICE

- Repeated autoclaving may cause the handpiece to become discolored due to heat. However, this is due to properties of the product and is not a problem in terms of quality.

2. Component Names



No.	Parts Name	Quantity
1	Control Unit	1
2	VA Bottle	2
3	AC Power Cord	1
4	Varios2 Handpiece (Optic or Non-Optic)	1
5	Handpiece Cord (Unshielded 2M) (Optic or Non-Optic)	1
6	Foot Control	1
7	Tip Wrench	3
8	Tip (G4, G6, G8)	1 each
9	O-Ring (Thin section) (For VA Bottle)	2
10	O-Ring (Thick section) (For VA Bottle)	2
11	O-Ring (For Handpiece Cord)	2
12	Water Connector (Option)	1
13	Water Tube (Option)	1
14	Spanner Wrench (5x8) (Option)	2
15	Tip Cover S (Option)	1
16	Tip-Power Guide	1*
17	Tip Card	1*
18	Operation Manual	1*

* These are not on photo above.


* Operation Principle

A sinusoidal electrical signal, at ultrasonic frequency ($f > 20\text{Khz}$), is delivered by the generator. This signal is applied to the 'piezoelectric ceramic' located inside the transducer. Piezoelectric ceramic converts this signal into mechanical vibrations. These vibrations are at the same ultrasonic frequency as the electrical signal. The mechanical vibrations are propagated towards the distal end of the transducer. The "TIP" insert, which is attached at the distal end of the transducer, vibrates at ultrasonic frequencies and makes it possible to achieve the aimed purpose.

3. Name and Function of each part

◆ Operation Panel and Display

* If you purchase the Optional products such as Water Tube and Water Connector you can use Tap Water.

Mode Select Keys

You can select Operation Mode to pressing this key. (Perio, Endo and General) The Control Unit can resume power level, water volume and irrigation mode for each Operation Mode.

Power Level Keys

Press key to select Power o/p Level . There are 11 levels (0 to 10). There is no output vibration at level 0 (zero). (Fig.1). The Bar Graph Indicator and Numerical Display will change simultaneously.


Fig.1

Irrigation Select Key

Press key to select 'R' or 'L' Bottle. Front panel display and Bottle Selection Indicator simultaneously change in position. Pressing the Irrigation Select Key for more than one second will select Tap Water Mode.

Auto-Cleaning Key

Press key to select Auto Cleaning Mode, For detail refer to 12. (5).

Bottle Water Adjustment Knob

Water Volume Adjustment can be made prior to the tip vibrating you can adjust the Water Volume during Bottle Irrigation or during a wait before Tip vibration start. If the setting is not applicable (too low or too high) for the Control Unit, it may beep. During the operation, Front Panel displays the current Power Level. However, keep turning the knob more than a second; it may change the Water Volume.

CAUTION

- Do not turn the knob fast. It may not sense the movement if it turn fast.
- Water Volume can set during 5ml/min to 45ml/min.
- Operation sound may different between Right and Left Bottle.
- During adjustment of water volume, Numerical Display indicate “—”.

Tap Water Adjustment Knob

You can adjust tap water volume by this Knob. (Even the tip vibration).

4. Prior to Operating System

(1) Water System Setup

• Use of Bottle

- 1) Remove the Dust Cover from the Bottle Base Connector. (Fig. 2)
- 2) Remove the cap of the VA Bottle and fill solution to the desired level.
- 3) Close the cap of VA Bottle and insert the Bottle Joint into the Bottle Base Connector until it clicks. (Fig. 3)

To remove the Bottle, pull it up.


Fig.2

CAUTION

- Use the VA Bottle Set 400 only for Varios 970.
- Before filling solution to the VA Bottle, check the Gasket inside the bottle cap is clean. (Fig. 4)
- Do not use a sharp tool to clean the Gasket or do not allow any impact on to the product. It may cause malfunction.
- Insert the Bottle straight. (Damage to seals may result).
- Keep the Gasket clean. When it becomes dirty by water or antiseptic solution, wash by clean water it immediately.
- The Gasket is consumable. *Optional Gasket: Order Code Z1047350


Fig.3

NOTICE

- The Bottle calibrations are printed on both sides of the Bottle and can be read accurately from the fill position or mounted on the Control Unit.
- Mount the Dust Cover when not in use.


Fig.4

• Use of Tap Water (Option)

- 1) Remove the Cover from the Tap Water Connector. (Fig. 5)
- 2) Connect the filter side of the Water Tube deep into the Control Unit Tap Water Connector (Fig. 6).
- 3) Connect the water tube to the water outlet on the Dental Unit.


Fig.5


Fig.6

CAUTION

Ensure the water from the dental unit runs clear before connecting to the scalar.

NOTICE

- Insert the Water Tube firmly into Control Unit.
- Pushing the White Ring, (the Quick Connector Release Ring) on the Tap Water Connector, gently pull the tube to remove. (Fig. 7)
- When the water tube is not connected, mount the cover on the Tap Water Connector.


Fig.7

(2) Foot Control Connection

Connect Foot Control Plug and the Control Unit with [▲] mark on the upper surface of the plug. (Fig.8)


Fig.8

(3) Handpiece Cord Connection

Insert Handpiece Cord Plug into Control Unit. [▲] Mark side is upper surface. Do not insert it up-side-down. (Fig.9)


Fig.9

CAUTION

Check that the Handpiece Cord Plug is clean & dry before connecting. (Fig.10)

Fig.10

(4) Mounting Power Cord

Insert the Power Cord into the AC Power Cord Connection Jack at the back of the Control Unit. (Fig. 11)

Fig.11

CAUTION

- Ensure Power is OFF on the Control Unit during the Power Cord Connection. It may cause Fuse to blow.
- Do not connect the cord in wall outlet before connecting AC Power Cord.
- Do not pull the AC Power Cord forcibly.
- Do not unplug the Power Cord or handpiece cord while pressing on the Foot Control.

5. Sterilization before first time use

NSK Varios components are not sold sterile and should be sterilized prior to first time you use. Autoclave sterilization is recommended.

Please note the following when sterilize components.

Follow your local directives, standards, and guidelines.

Sterilize components in a U.S. FDA-cleared sterilization pouch, complaint to ISO 11607-1.

Ensure that you only remove dry sterile goods from the sterilizer.

Store sterile goods dust-free and dry.

NOTICE

- NSK recommends sterilization according to ANSI/AAMI ST79 and EN 13060 (class B).
- Follow the manufacturer's instructions for use that came with your sterilizer.
- For more information about sterilization standards, refer the medical device regulations specific to your location.

■ Autoclave Procedure

- 1) Pick a tip up from the tip case. Take handpiece out of the packing bag before sterilization.
- 2) Insert sterilizable components into a U.S. FDA-cleared sterilization pouch, complaint to ISO 11607-1. Seal the pouch.
- 3) Autoclave conditions are below:

< When sterilizing a Tip>

Steam Sterilization Cycle					
Type	Gravity Displacement		Pre-Vacuum (Dynamic Air Removal)		
Minimum Temperature	132°C	135°C	132°C	134°C*	135°C
Full Cycle Time	15 minutes	10 minutes	4 minutes	3 minutes	3 minutes
Minimum Drying Time	15 minutes	30 minutes	20 minutes	30 minutes	16 minutes

* Not for use in the U.S.

<When sterilizing components except for a Tip>

Steam Sterilization Cycle				
Type	Gravity Displacement			Pre-vacuum
Minimum Temperature	121°C	132°C	135°C	135°C
Full Cycle Time	20 minutes	15 minutes	10 minutes	10 minutes
Minimum Dry Time	10 minutes	10 minutes	10 minutes	15 minutes

CAUTION

- Do not exceed 138°C during sterilization.

- 4) Keep the components in the autoclave pouch until ready to use.
- 5) Store in a dry and dust-free location up to the period specified by the sterilization pouch manufacturer. If sterility cannot be confirmed, sterilize again prior to use.

6. Mounting and Removing the Handpiece

Align the Dots on the Handpiece and the Handpiece Cord. Push handpiece into connector.

To remove the handpiece, grip the Handpiece and Handpiece Cord and pull to part handpiece and cord. (Fig. 12)

WARNING

To avoid Electrical Shock Do not touch the handpiece backend electrical contacts.


Fig.12

CAUTION

- Before connecting the handpiece cord to the handpiece, face the handpiece connector side so light shines into it and check the end of the handpiece. If the light transmitted by the Optic Fiber End Face at the end is uniformly bright, it is functioning normally. If part or all is dark, the Optic Fiber End Face may be cracked or broken.
- Always confirm that the handpiece is correctly seated and locked into place.
- Do not connect or use Handpiece other than included one (Varios2 handpiece).

7. Mounting and Removing Tip

- 1) Turn Tip lightly by hand, and install it. (Fig. 13)
- 2) Tip will insert from the bottom hole of Tip Wrench. Align the four corner of the Tip base area into the four corner of Tip Wrench. And turn it clockwise until it clicks. (Fig. 14)

* Do not touch the top part of Tip to avoid an injury. (There is the case that is longer than height of TIP WRENCH)

To remove the Tip, turn it counterclockwise with the Tip Wrench.


Fig.13


Fig.14

Caution for Tip Usage

- Check the Tip before use. (Flush, Damage, Bending or Rust)
- Do not exceed Maximum Power Level for Tip. Damage to tooth structure and Tip may result.
- Do not hit ceramic prosthesis with Tip during scaling. Tip Damage may result.
- Do not hit metal or prosthetic crown except for removing them. Tip could break and fall into mouth.
- Do not hit gingival, mucosa and/or skin. It could cause damage and/or burn injury.
- During its operation, the whole tip vibrates due to ultrasonic waves. Do not bring any part of the tip into contact with soft tissue, gums or skin. Regardless of use of coolant water, the tip in vibration could cause burn injury.
- The “dry tip*” (a tip which does not require coolant water) utilizes heating to some extent for effective treatment. During or after its operation, the whole of such tip could be overheating. Do not bring any part of the tip into contact with soft tissue, gums or skin. Ensure that the power setting of the Control Unit does not exceed the maximum power written on the Tip Case label. Adjust the power from the lowest level. During use of the tip, ensure that burn injury does not occur in or around the treatment area. Also ensure that the tip does not affect the treatment area and oral cavity.
* E5, E6, E7, E8, G21, G22, G28
- Do not sharpen and/or bend the Tip. Tip may damage and not generate enough vibration during scaling.
- During cutting, Tip will gradually wear away, as the Tip wears the stroke will get smaller and decrease cutting efficiency. When level drops too far, change the Tip. (tip card check)
- DO ENSURE When securing tip to use the tip wrench as supplied, inefficient cutting will result.
- DO ENSURE before attaching Tip, Cleanliness of the tip screw, inefficient cutting will result.
- To avoid personal injury DO ENSURE tip is removed prior to disconnecting the handpiece from the handpiece cord.
- If you feel the Tip is not vibrating, remove it from an operative site, and press the Foot Control again. If this does not improve the condition, Ensure the tip is secure, turn the power off and restart it.
- When mounting the Tip, always use grooves and Tip Wrench as supplied.
- Ensure that water volume must be “0”, when you use Tip which does not appear of water.
- Tip Wrench is consumable. For reliable operation replace annually.

8. Operating Procedures

(1) Water System Setup

CAUTION

Do not spray using ozonated water, high oxidization potential water (strongly acidic water, superacidic water, electrolyzed water, etc.), or hypochlorous acid water. Doing so may lead to clogging, corrosion of metal parts, etc.

• Use of Bottle

- 1) Check that the VA Bottle is filled to the proper level.
- 2) Make sure that the cap of the Bottle is secure and not leaking.

CAUTION

DO ENSURE liquid temperature is below 35°C.

• Use of Tap Water

- 1) Ensure water tube is firmly connected.
- 2) Open the dental unit's water valve. (Set water pressure between 0.1-0.5MPa (1-5 kgf/cm²)).

CAUTION

Ensure that the water purity is adequate for the specific, intended procedure.

(2) Power On

Connect the AC Power Cord to the wall outlet. Turn the Power Switch on the Control Unit. Front Display will illuminate.

Fig.15

(3) Power Level Setting

DO ENSURE Power setting does not exceed the recommended Power Level (Tip-Power Guide included in the package.)

1) Select the Operating Mode with the Mode Select Keys on the Front Panel. The Indicator light over the selected mode will illuminate. (Fig. 16)



Fig.16

* Power Level for each mode

2) Set the power level with the Power Level Key on the Front Panel. The Bar Graph Indicator and Numerical Display will indicate the selected power level. (Fig. 17) Make sure the power level is set in the appropriate range for the attached Tip.

NOTICE

- Press & Hold the Power Level Key will increase or decrease the Power Level.
- If the Power Level is 0 (zero) and set the water volume, Tip will not oscillate but water comes out from the handpiece.


Fig.17

(4) Irrigation setting

Select the Irrigation Mode (L Bottle, R Bottle or Tap Water) with the Irrigation Select Key on the Front Panel. (Fig. 18) The Indicator light over the selected mode will illuminate. Press & Hold the Irrigation Select Key to select Tap Water Mode.


Fig.18

(5) Operate Varios 970 / 970 LUX

By stepping on the Foot Control, the tip vibrates and spraying starts (except for tips that do not spray) and the handpiece LED lamp lights up (Varios 970 LUX). When you remove your foot from the foot control, tip vibration and water spraying stop and the handpiece LED lamp turns off. (Varios 970 LUX).

• Water Supply Volume Adjustment

Turn the Water Adjustment Knob counterclockwise gradually to increase the supply volume. (Fig. 19) For detail, refer to P6 Bottle Water Adjustment Knob or P7 Tap Water Adjustment Knob.

Fig.19

CAUTION

- While pressing the Foot Control and switching the power 'ON'. The Control Unit will display "F" and sound a beep, for your safety the Control Unit will not operate. Remove your foot from Foot Control to cancel and reset.
- Bar Graph Display (Fig.20)
 - Minimum Irrigation -> One white and blue LED.
 - No Irrigation -> Blue LED only
- Always use the water supply. If water supply is insufficient, handpiece will overheat and patient's tooth surface can be injured.
- Verify that the water spray is clean and of adequate volume before use.
- If irrigation volume set low, sometimes irrigation water is difficult to come out from the Tip. When it happened, set volume again after setting up high volume.
- During Water Adjustment Knob operation;
 - Numerical Display: Display "-"
 - Bar Graph: Display current volume of water


Fig.20

(6) After the treatment

Release the Foot Control and Power off the Control Unit.

• Use of Bottle

Thoroughly wash the Bottle (s) Water Supply system. Refer to 12. (5) Auto Cleaning (Cleaning of Irrigation Tube).

CAUTION

When using medicated solutions, clean the entire Irrigation System thoroughly.

• Use of Tap Water

Close the dental unit's water valve.

NOTICE

- LED of the handpiece will remain 'On' for approx 5 seconds after Foot Control is released. (Varios 970 LUX)
- When the Control Unit is Power off, the last mode settings in use are automatically retained in memory.

◆ Initialized Program (Factory Setting)

Press both Auto-Cleaning and Power Key on the Control Unit to initialize the Factory memory Setting. Do not release Auto-Cleaning Key until the beeping sound from the Control Unit. (Initial Mode is Perio)

	Power	Flow amount (L, R each)	Irrigation Mode	Initial Mode
Perio	1	10	L Bottle	●
Endo	1	10	L Bottle	
General	1	10	L Bottle	

◆ During the Handpiece operation :

Possible: Power Level and Water Volume adjustment.

Impossible: Operation Mode and Irrigation Mode setting, Auto Cleaning.

9. Provided Scaler Tips

G4

The end of the Tip is thin and for supragingival fine scaling and interdental scaling. The round cross-section allows tooth surfaces to be finished without causing damage.

Apply the top of the Tip on the tooth plane and move it sideways finely in the same way as G8 Tip. (Fig. 21)

Fig.21

G6

Removal of supra and subgingival calculus. It provides easy access to interdental spaces and narrow pockets.

Insert the top of the Tip into the periodontal pocket and move it slowly. The top of the Tip is sharp so that it could remove tartar on long coroner and retracted gingival. (Fig. 22)
Clean periodontal pocket at low power.

Fig.22

G8

Removal of supragingival and interdental calculus. This Tip can be used in all quadrants and is very useful for the removal of hard calculus.

Apply the top of the Tip on the tooth plane and move it sideways finely along the neck of tooth. (Fig. 23)

Fig.23


CAUTION

Tip is article of consumption. We recommend periodical replacement. About time of replacement, check the Tip Card.

◆ How to use the Tip Card

- 1) Place the neck of the Tip in the cut out.
- 2) Check wear of the Tip.
- 3) See the green, yellow and red line to check wear of the Tip. *See below what each color means. At NSK we recommend to replace a Tip when the Tip meets the yellow line (wear of 1mm) to guarantee safe and effective use.

Tip Card

Green: No wear - Tip is OK

Tip replacement is not necessary.

Yellow: Wear of 1mm - Tip is showing some wear

Tip replacement is recommended.

Red: Wear of 2mm - Tip is badly worn

Tip replacement is necessary.

! CAUTION

Tips are consumables. The efficiency of dental scaling decreases approximately 25% when the top of the Tip wears 1 mm and approximately 50% when it wears 2 mm. In addition, the vibration condition changes owing to the wear, which may damage a patient's tooth surface. Check the Tip wear condition with the Tip Card periodically, and replace the Tip with a new one in good time.


Fig.24

10. How to Use Tip Cover S (Option)

Grip the Tip Cover S and insert it to the Tip.

To remove, grip the Tip Cover S and the handpiece & pull. (Fig. 25)

! CAUTION

- Carefully insert the Tip into the Tip Cover S. Avoid injuring the fingers.
- The Tip Cover S is not designed for use as a Tip changing tool.



Fig.25

11. Handpiece Holder

While the Handpiece is not in use, put the Handpiece in the Handpiece Holder.

The Handpiece Holder is adjustable. (Fig. 26)

CAUTION

Do not put excessive load to the Handpiece Holder to prevent from breaking down and deformation.

NOTICE

To prevent injury, always mount Scaler Tip Cover (S).

Fig.26

12. Care and Maintenance

(1) Cleaning of Optic Fiber (Varios 970 LUX)

Wipe the debris off the end of the Optic Fibers at the handpiece with alcohol soaked cotton swab. (Fig. 27)

CAUTION

Do not use any sharp pointed tools to clean the Optic Fiber End Face. In case the light degradation, contact your dealer.


Fig.27

(2) Cleaning the Handpiece cord

Remove the handpiece after use on each patient and clean it as described below.

- 1) Wipe the surface of the handpiece cord with a cloth soaked in alcohol.
- 2) Carefully wipe the handpiece cord plug with an alcohol-immersed cotton swab. If it is difficult to use a cotton swab, carefully wipe with a towelette wound around a thin stick-shaped object.


Fig.28

CAUTION

Do not use a sharp pointed stick or push the terminal part when cleaning the handpiece cord plug. Doing so may cause damage, resulting in contact failure (Fig. 28).

(3) Changing O-Ring

- **Handpiece Cord**

An O-Ring is located in the Handpiece Cord Connector. Use a pointed tool to remove, and mount new O-Ring into the groove. (Fig. 29)

* Optional O-Ring: Order Code D0310020080

Fig.29

- **VA Bottle**

Remove two O-Rings at the Bottle Joint with a pointed tool, and mount new O-Rings into the grooves. (Fig. 30)

* O-Ring (Thick section) : Order code D0310075150

O-Ring (Thin section) : Order code D0312090100


Fig.30

(4) Changing the Irrigation Pump

- 1) Remove the Bottle, the Power Cord, the Handpiece Cord and the Foot Control from the Control Unit.
- 2) Turn back the Control Unit. Hang a finger on "A" point and pull up the bottom cover to remove.

Fig.31

Picture below is shows inside of the Control Unit.




Fig.32

- 3) Remove the Irrigation Tube from the Control Unit. (Bottle side and Front Panel side.) (Fig. 33, 34)
- 4) Remove the Connector Ring from the Irrigation Tube. Do not dispose it. You can use the Rings to the replacement Irrigation Pump.
- 5) Turn the Irrigation Pump counterclockwise until it clicks and pull it out. (Fig. 35)
- 6) Mount the Connector Ring to the new Irrigation Pump. Observing Ring direction. (Fig. 36)
- 7) Align the replacement Irrigation Pump with the Drive Shaft. Turn clockwise until it clicks. (Fig. 35)
- 8) Mount the Irrigation Tube opposite procedure of removing (Fig.33). Connector Ring should firmly into the Control Unit until it stops. (Fig.37)

* Bottle side

Fig.33

* Front Panel side

Fig.34


Fig.35



Fig.36

* Front Panel side

Fig.37

- 9) Align the Bottom Cover Hook and hole on the Control Unit. Mount the Bottom Cover.

* Optional Irrigation Pump: Order Code 10000643 (Not included the Connector Ring.)

CAUTION

- If water is spilled out the irrigation pump, wipe it off and allow drying completely prior to use. If water gets inside the irrigation pump the roller may slip and fail to pump.
- Before replacing the Irrigation Pump, wipe off excess water on pump and Drive Shaft. The wet drive shaft and rollers can be slippery and cause improper operation.
- Wipe dirt and water off the Drive Shaft from bottom up. (Fig.38)
- Insert the replacement Irrigation Pump into the Drive Shaft straight (slow and soft) to prevent damaging rollers in pump.
- Run the replaced Irrigation Pump about 10 seconds on largest setting of Water Volume before operation to adopt Irrigation Tube to new pump.
- Ensure Irrigation tube has no kink or twists If tube is set incorrectly, Irrigation Water may not come out.
- Do not pull the tube when the bottom cover is closed.


Fig.38

NOTICE

- Perform periodical cleaning for the Drive Shaft with socked alcohol cloth. Dirt on Drive Shaft may cause an incorrect pump operation.
- The pump is consumable. If the irrigation volume decreases markedly, replace pump.

(5) Auto Cleaning (Cleaning of Irrigation Tube (Use of Bottle))

NOTICE

- The Auto Cleaning function is not a device disinfection step. The Auto Cleaning function is intended only to purge and rinse the system of any remaining solution used during a dental procedure.
- After each use, remove all of the disinfectant solution and perform the "Auto Cleaning" (or manual cleaning) procedure. If the system is not cleaned, residual disinfectant can cause damage to the system.
- During Auto Cleaning, water comes out from the handpiece. Perform cleaning after turning handpiece into a cup.

- 1) Remove the 2 Bottles from the Control Unit.
- 2) Clean inside of the Bottle.
- 3) Half fill the bottle with purified water (DO NOT USE SALINE)

CAUTION

Use only distilled water for cleaning.

- 4) Install the cap on the Bottle. Install the Bottle Joint into the Bottle Base Connector. Install it until clicks into place. Improper connection may cause water leakage. Make sure the connection is tight.

CAUTION

- Perform Auto-Cleaning without tip.
- Make sure the handpiece and handpiece cord are firmly attached.


- 5) To perform the Auto-Cleaning, keep pressing the Auto-Cleaning Key (more than 1 second). It takes 30 seconds per bottle to clean. The Numerical Display will alternately display "A" and "C", the Bar Graph displays time remaining. Single display (Bar Graph Display) is 6 seconds. When five displays of Bar Graph disappeared, Bottle will be changing to the other side.

To cancel the Auto-Cleaning, press Auto-Cleaning Key twice.

NOTICE

During Auto-Cleaning, LED of the handpiece does not illuminate. (Varios 970 LUX)

- 6) When the Auto-Cleaning is finished, the Control Unit returns to the settings prior to cleaning. Remove the both bottles from the Control Unit by pulling straight up. Clean thoroughly, rinse and dry.

◆ Following method is also available for cleaning. (Manual Cleaning)

- 1) Remove the Bottle from the Control Unit.
- 2) Open the cap of the cleaned Bottle and fill it with distilled water.
- 3) Close the cap firmly and insert the Bottle Joint into the Bottle Base Connector on the Control Unit until it clicks.
- 4) Operate the Control Unit about 30 seconds with water supply at maximum setting.

NOTICE

The Control Unit does not perform in Auto-Cleaning in Tap Water.

(6) Changing Water Filter (Option)

If you use Tap Water, change the Water Filter as it may necessary.

- 1) Close the water valve of the dental unit.
- 2) Mount two Spanner Wrenches (5x8) and turn those as shown in Fig.40.
- 3) When the Water Filter case is separated, the Water Filter can be removed as shown in Fig. 41.
- 4) Replace with new (Order Code U387042) and reassemble the filter in the reverse order.


Fig.40


Fig.41

13. Periodical Maintenance Checks

In order to use the product properly and safely, perform operation checks once a year following the steps shown in "4. Prior to Operating System", "6. Mounting and Removing the Handpiece", "7. Mounting and Removing Tip", and "8. Operating Procedures". Contact your Authorized NSK Dealer if any abnormalities are found.

14. Sterilization

- All sterilizable components of this device must be sterilized after each patient.
- The recommended method of sterilization is by autoclave sterilization.
- ONLY the Tip (except diamond coated scaler tip), Handpiece, Tip Wrench, Tip Holder and Tip Cover S can be autoclaved. The diamond coated Scaler Tip is to be sterilized prior to first use, then disposed of as this Tip is single use only.

NOTICE

- NSK recommends sterilization according to ANSI/AAMI ST79 and EN 13060 (class B).
- Follow the manufacturer's instructions for use that came with your sterilizer.
- For more information about sterilization standards, refer the medical device regulations specific to your location.

 The surface of the handpiece can be washed and disinfected with a Thermo-Disinfector prior to Sterilization.

■ Autoclave Procedure

- 1) Remove the Tip after use. (Refer to 7. Mounting and Removing Tip)
- 2) Clean all the surfaces of the components using a small tooth brush with ethanol, or ethanol with isopropyl alcohol to remove dirt and debris then dry thoroughly.
- 3) Insert sterilizable components into a U.S. FDA-cleared sterilization pouch, compliant to ISO 11607-1. Seal the pouch.
- 4) Autoclave conditions are below:

< When sterilizing a Tip >

Steam Sterilization Cycle					
Type	Gravity Displacement		Pre-Vacuum (Dynamic Air Removal)		
Minimum Temperature	132°C	135°C	132°C	134°C*	135°C
Full Cycle Time	15 minutes	10 minutes	4 minutes	3 minutes	3 minutes
Minimum Drying Time	15 minutes	30 minutes	20 minutes	30 minutes	16 minutes

* Not for use in the U.S.

<When sterilizing components except for a Tip >

Steam Sterilization Cycle				
Type	Gravity Displacement			Pre-vacuum
Minimum Temperature	121°C	132°C	135°C	135°C
Full Cycle Time	20 minutes	15 minutes	10 minutes	10 minutes
Minimum Dry Time	10 minutes	10 minutes	10 minutes	15 minutes

CAUTION

- Do not exceed 138°C during sterilization.

- 5) Keep the components in the autoclave pouch until ready to use.
- 6) Store in a dry and dust-free location up to the period specified by the sterilization pouch manufacturer. If sterility cannot be confirmed, sterilize again prior to use.
- * The recommended sterilization procedures require the use of FDA-cleared sterilizers, sterilization trays, sterilization wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended.

CAUTION

- Do not sterilize by ultraviolet ray. The handpiece could discolor.
- If autoclaved with other instruments stained with chemical solution, it could strip the plating and make the surface black.
- Do not autoclave any parts (the Control Unit, AC Power Cord, Foot Control, Handpiece Cord, Tip Card, Tip Power Guide, O-Ring). Other than those that can be subjected to autoclave sterilization. Perform alcohol disinfection to the Control Unit, AC Power Cord, Foot Control, Handpiece Cord, Tip Card, Tip-Power Guide including after every patient.
- Do not wipe with, or clean or immerse in, high acid water or sterilizing solutions.
- Diamond coated Scaler Tip is for single use only. Do not sterilize and use again.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Do not touch the product immediately after autoclaving may cause injury.

NOTICE

Repeated autoclaving may cause the handpiece to become discolored due to heat. However, this is due to properties of the product and is not a problem in terms of quality.

15. Periodical Maintenance Checks

Refer to the check sheet below and perform periodical maintenance checks every one year. Contact your Authorized NSK Dealer if any abnormalities are found.

Points to Check	Details
Handpiece	Operate the Handpiece to confirm there is no abnormal vibration, noise, heat, etc. Also check that the LED lights (Varios 970 LUX) and water flows.

16. Troubleshooting

When trouble is found, please check the followings prior to consulting your dealer.

Problem	Probable Cause	Cause	Solution
No / Poor vibration.	The Front Panel does not light, even if the Power Switch is ON.	The Power Cord or the Jack is disconnected.	Correctly insert the Power Cord or the Jack.
		The Fuse is burned out.	Contact dealer.*
	The Tip does not generate vibration, in spite of depressing the Foot Control.	The Tip is not tightened firmly.	Tighten the Tip until the Tip Wrench clicks.
		Worn Tip.	Replace the Tip.
		Power has not been correctly adjusted for the Tip.	Adjust the power on the Power Guide or Tip case label. Do not exceed.
		The Foot Control is disconnected.	Connect the Foot Control correctly.
		Failure of vibrator in the handpiece.	Contact dealer.*
Failure of internal components of the Foot Control.	Contact dealer.*		
The Tip is bent or broken.	—	Power has not been properly adjusted for the Tip.	Adjust the power level the Power Guide or Tip case label. Do not exceed.
The Tip is flying away.	—	The Tip is not tightened firmly.	Tighten the Tip until the Tip Wrench clicks.
Noise from the handpiece.	—	Power has not been properly adjusted for the Tip.	Adjust the power level on the Power Guide or Tip case label. Do not exceed.
		The Tip is not tightened firmly.	Tighten the Tip until the Tip Wrench clicks.
		Failure of vibration in the handpiece or the Control Unit.	Contact dealer.*
The handpiece is overheating.	—	Power has not been properly adjusted for the Tip.	Adjust the power level on the Power Guide or Tip case label. Do not exceed.
		The Tip is not tightened firmly.	Tighten the Tip until the Tip Wrench clicks.
		Failure of vibration in the handpiece or the Control Unit.	Contact dealer.*
No Irrigation supply and/or unstable Irrigation supply (Use of Bottle)	The Irrigation Pump is running.	The tube twisted.	Straighten the twisted Irrigation Tube.
	The Irrigation Pump is stopping.	Time to replace Irrigation Pump. (Approx. 500hours after used.)	Replace with new Irrigation Pump (Refer to 12. (4) Changing the Irrigation Pump).
No / Poor water. (Use of Tap Water)	The water does not reach to the Control Unit.	—	Check the water circuitry and supply to the Control Unit. Water pressure : 0.1-0.5MPa (1-5kgf/cm ²)
		The Water Adjustment Knob is closed.	Turn the Water Adjustment Knob and adjust to the appropriate volume.
	Check to see if water reaches the Control Unit.	Disconnected Irrigation supply at low volume range. (less than 10ml/ min.)	No problem. Turn the Water Adjustment Knob and increase the Irrigation volume.
		The Water Filter is clogged.	Replace with new Water Filter (Refer to 12. (6) Changing Water Filter (Option)).

Problem	Probable Cause	Cause	Solution
Water leakage.	Water is leaking from the joint between the Irrigation Tube and the Irrigation Connector.	The Irrigation Tube is not connected correctly.	Firmly insert the Irrigation Tube into the Irrigation Connector inmost.
	Water is leaking from the joint between the handpiece and the cord.	O-Ring at the handpiece cord is worn or damaged.	Replace with new O-Ring (Refer to 12. (3) Changing O-Ring •Handpiece Cord).
	Water is leaking from the Control Unit.	The water circuitry in the Control Unit is damaged.	Contact dealer.*
Handpiece LED does not illuminate. (Varios 970 LUX)	Tip oscillates, but Handpiece LED turns on and off.	The handpiece is not connected into the Handpiece Cord correctly.	Firmly insert the handpiece into the Handpiece Cord inmost.
	Tip oscillates, but Handpiece LED does not turn on.	Disconnection in the Handpiece Cord, or failure in the Control Unit.	Contact dealer.*
Start Beeping	Beeping while power on.	Depress Foot Control.	Release the Foot Control.
	Beeping while stopping vibration of Tips.	Abnormal heating of the Control Unit.	Stop the operation and leave until Control Unit becomes cool.

* Repairs cannot be made by the customer.

17. Protection Circuit

It may overheat inside when you use this Control Unit in more than Power 8 at G mode for long time.

In this case, Protection Circuit reduces the Power automatically. (Power 7)

Bar Graph Indicator from 8 to 10 flashes. (Fig.42)

After Protection Circuit is released, the flashes stop.

However, Power Level can not automatically increase. If needed, increase manually.


Fig.42

NOTICE

- During Protection Circuit function (during Bar Graph Indicator flash), the Control Unit can not increase the Power Level.
- If Power Level decreases less than 7, Bar Graph Indicator stops flashing. However, it the Power increase more than 8, flashes it again.

18. Error Code

If an operational problem occurs numerical Display shows the error code to allow an immediate problem diagnosis.

Error Code	Error	Check / Remedy
E 0	Self-Check Error	Contact dealer.
E 1	Circuit Failure	Contact dealer.
E 7	Does not vibrate	Contact dealer.
E 9	Handpiece Self Check Error	Confirm connection of the handpiece. Power on the Control Unit again. Leave the Control Unit until it become cool down and powers it again. When an error can not be eliminated, Contact dealer.
E 10	Circuit Failure	Contact dealer.

**"E" and the number alternately display on the Display.

19. Specification

Type	NE255
Power Source	AC120V 50/60Hz AC230V 50/60Hz
Vibration Frequency	28-32kHz
Maximum Output	11W
Rated Power	29VA
Water Pressure	0.1-0.5MPa (1-5kgf/cm ²)
Lighting	Varios 970 : No Varios 970 LUX : Yes
Bottle Volume	400mL (Per Bottle)
Dimensions	W160 x D270 x H190mm (Including Bottle)
Weight	2.1kg (Except attachment)
Use Environment	Temperature 0 - 40 °C (The liquid must not freeze up) Humidity 30 - 75 % Atmospheric pressure 700 - 1060 hPa
Store Environment	Temperature -10 - 60 °C Humidity 10 - 85 % Atmospheric pressure 500 - 1060 hPa

20. Spare Parts

Model	Products	Oeder Code	Model	Products	Oeder Code
VA Bottle Set 400		Z1047002	VA Bottle Cap 400		10000652
VA Bottle 400		20000947	Gasket		Z1047350
Water Tube Set		U387040	Irrigation Pump		10000643
Water Connector		U387030	Tip Wrench (CR-10)	 	Z221076
Water Filter		U387042	Tip Holder	 	Z221A080
Spanner Wrench (5x8)		Y1001301	Tip Cover S	 	Z217851
O-Ring (for Handpiece Cord)		D0310020080	O-Ring (Thick section) (For VA Bottle)		D0310075150
			O-Ring (Thin section) (For VA Bottle)		D0312090100

 Autoclavable at 135°C max.

21. Disposing product

Consult with dealer from whom you purchased it about waste disposal.

22. Warranty

Manufacturer warrants its products to the original purchaser against defects in material and workmanship under normal practices of installation, use and servicing. Such expendable items as O-Rings and Irrigation Pump are not covered by this warranty.

Symbols

TUV Rhineland of North America is a Nationally Recognized Testing Laboratory (NRTL) in the United States and is accredited by the Standards Council of Canada to certify electro-medical products with Canadian National Standards.

Follow the waste of electric and electronic equipment (WEEE) Directive (2002/96/EC) to dispose of the product and accessories.

Consult operation instructions.

Manufacturer.

This conforms to CE European Directive of "Medical equipment directive 93/42/EEC."

Type BF applied part.

Authorised representative in the European community.

Protected against vertically falling water drops.

Autoclavable up to Max.135°C. *for detail see Sterilization.

The surface of the product can be washed and disinfected with a Thermo-Disinfector prior to Sterilization.

Marking on the outside of Equipment or Equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment.

Single use. Do not reuse.

Rx Only Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

GS1 DataMatrix for Unique Device Identifier.

Guidance and manufacturer's declaration - electromagnetic emissions		
The Varios 970 / Varios 970 LUX is intended for use in the electromagnetic environment specified below. The customer or the user of the Varios 970 / Varios 970 LUX should assure that is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR11/EN55011	Group 1 class B	The Varios 970 / Varios 970 LUX uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Varios 970 / Varios 970 LUX is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supply network that supplies buildings used for domestic purposes.
RF emissions CISPR11/EN55011	Group 1 class B	
Harmonic emissions EN/IEC61000-3-2	class A	
Voltage fluctuations/flicker emissions EN/IEC61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The Varios 970 / Varios 970 LUX is intended for use in the electromagnetic environment specified below. The customer or the user of the Varios 970 / Varios 970 LUX should assure that it is used in such an environment.			
Immunity test	EN/IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) EN/IEC61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst EN/IEC61000-4-4	±2kV for power supply lines ±1kV for input/output	±2kV for power supply lines ±1kV for input/output	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC61000-4-5	±1kV line(s) to line(s) ±2kV line(s) to earth	±1kV line(s) to line(s) ±2kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN/IEC61000-4-11	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 secs	<5% Ut(>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% Ut (>95% dip in Ut) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Varios 970 / Varios 970 LUX requires continued operation during power mains interruptions, it is recommended that the Varios 970 / Varios 970 LUX be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field EN/IEC61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: Ut is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic immunity			
The Varios 970 / Varios 970 LUX is intended for use in the electromagnetic environment specified below. The customer or the user of the Varios 970 / Varios 970 LUX should assure that it is used in such an environment.			
Immunity test	EN/IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF EN/IEC61000-4-6	3Vrms 150 kHz to 80MHz	3Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Varios 970 / Varios 970 LUX, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80MHz to 800MHz $d = 2.3\sqrt{P}$ 800MHz to 2.5GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF EN/IEC61000-4-3	3V/m 80MHz to 2.5 GHz	3V/m	
NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobiles radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Varios 970 / Varios 970 LUX is used exceeds the applicable RF compliance level above, the Varios 970 / Varios 970 LUX should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Varios 970 / Varios 970 LUX.			
b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3 V/m.			

Cables and accessories	Maximum length	Complies with	
Handpiece cord	2 m	RF emissions, CISPR11, EN55011	Class B/ Group 1
Foot Control	2.5 m	Harmonic emissions, Voltage fluctuations/ flicker emission, Electrostatic discharge (ESD) Electric fast transient / burst Surge Voltage dips, short interruptions and voltage variations on power supply input lines Power frequency(50/60Hz) magnetic field Conducted RF Radiated RF	EN/IEC61000-3-2 EN/IEC61000-3-3 EN/IEC61000-4-2 EN/IEC61000-4-4 EN/IEC61000-4-5 EN/IEC61000-4-11 EN/IEC61000-4-8 EN/IEC61000-4-6 EN/IEC61000-4-3

Recommended separation distances between portable and mobile RF communications equipment and the Varios 970 / Varios 970 LUX.			
The Varios 970 / Varios 970 LUX is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Varios 970 / Varios 970 LUX can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Varios 970 / Varios 970 LUX as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d=1.2\sqrt{P}$	80MHz to 800MHz $d=1.2\sqrt{P}$	800MHz to 2.5GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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2019.02.04 001 ©
NSK America Corp www.nskdental.com