Multi Pad2 NLZ Built-in Motor System

INSTALLATION GUIDE

Rx Only

English1

Français.....37

TABLE OF CONTENTS

English

1	Preface.....	2
1.1	User and indications for use.....	2
1.2	Introduction.....	2
1.3	Symbol.....	3
1.4	Disposing of products.....	4
2	Package Contents.....	5
3	Specifications	7
3.1	Multi Pad2.....	7
3.2	NLZ BF.....	8
3.2.1	NLZ BF Control Board.....	8
3.2.2	AC – DC Adapter Board	9
3.2.3	Selector board.....	10
4	Mounting Multi Pad2.....	11
4.1	Mounting MP2 STAY	11
4.2	Mounting Angle Plate	13
4.3	Mounting Ferrite Core	14
5	Connection	15
5.1	Multi Pad2 and NLZ BF	15
5.1.1	Connecting Harness Assy to NLZ BF.....	15
5.1.2	DIP switch setting	16
5.1.3	NLZ BF speed control settings.....	17
5.2	Multi Pad2.....	19
6	Power Supply Connection	21
7	Cautions When Installing to Delivery Table.....	22
8	EMC Information (Electromagnetic Compatibility Information).....	23
9	Wiring	26
9.1	Wiring Diagram	26
9.2	Connecting joint	36

1 Preface

1.1 User and indications for use

- User
Dentist, Dental hygienist
- Indications for use
The NLZ Built-in Motor System is indicated for use in the field of prophylaxis dentistry, restorative applications including cavity preparation and endodontic therapy, prosthodontics applications such as crown preparations.
- Target for this guide
Manufacturers, distributors, and dealers of a dental chair table

1.2 Introduction

The NLZ Built-in Motor System is featured to integrate with a micromotor (NLZ BF) in one system. This guide is for installation and setup of the Multi Pad2 and micromotor (NLZ BF).

CAUTION

- Avoid possible damage to the product before installation as follows:
 - Disconnect power
 - Discharge electric static
 - Make sure the connection of the pin connector is correct
- If the Multi Pad2 is not activated or the power turns ON and OFF, check for a potential problem in the connection of cable connectors.
- When installing the motor cord, make sure the cord is safely and securely tied to the dental chair table. When securing the cord, make sure the inside tubes are not squeezed. If the cord is secured too tight, the air and water system will not function properly. Use the metal portion of the cord when installing the motor cord.
- Do not bend the water tube and air tube.
- This system complies with IEC 60601-1 and IEC 60601-1-2. Also confirm compliance with IEC 60601-1 and IEC 60601-1-2 in the complete product state with this system built-in.
- There is a risk of malfunction in an electromagnetic environment that exceeds the IEC 60601-1-2 test level.

1.3 Symbol

Consult operation instructions

Caution

Type B applied part

Marking on the outside of Equipment or Equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment

IPX1

Protected against vertically falling water drops

Dispose of this device and its accessories via methods approved for electronic device and in compliance with the Directive 2012/19/EU

Conforms to CE European Directive of “Medical device directive 93/42/EEC”

TUV Rheinland of North America is a Nationally Recognized Testing Laboratory (NRTL) in the United States and is accredited by the Standards Council of Canada to certify electro-medical products with Canadian National Standards

Distributor

Manufacturer

Medical Device

Caution: U.S. Federal law restricts this device to be used by or on the order of a licensed dental professional

Device Name for Health Canada

Catalog number (Order Code)

Serial No.

Date of manufacture

Unique Device Identifier

GS1 DataMatrix for Unique Device Identifier

Temperature limitation

Humidity limitation

Atmospheric pressure limitation

Keep dry

Fragile / handle with care

This is the correct upright position of the distribution packages for transport and/or storage


1.4 Disposing of products

In order to avoid the health risks to operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist must confirm the equipment is sterile.

Ask specialist firms who are licensed to dispose of specially controlled medical wastes to dispose the product for you.


2 Package Contents

NLZ BF


No.	Part Name	Model	Order Code
(1)	NLZ BF	NE348	U1104002
(2)	Motor	NLZ	E1152051
(3)	Motor Cord	NLZ CDB	EA1152X061
(4)	Dimension Conversion Joint Set	JOINT SET	-

Multi Pad2 / Angle Plate

No.	Part Name	Model	Order Code
(1)	Multi Pad2	NE344	Y1500471
(2)	Angle Plate	-	-
(3)	Binding Head Screws M2.6-L6	-	-
(4)	NLZ POWER HARNESS	-	-
(5)	Ferrite Core	-	-
(6)	UNITY BAND	-	-

MP2 STAY


No.	Part Name	Model	Order Code
(1)	Mounting stay	-	-
(2)	Clamp holder	-	-
(3)	Post A	-	-
(4)	Post B	-	-
(5)	Adjustment screw	-	-
(6)	Fixing nut	-	-
(7)	Mounting screws & nuts	-	-
(8)	Screw with washer	-	-
(9)	Cover	-	-
(10)	Cable tie	-	-

3 Specifications

3.1 Multi Pad2

Name	Multi Pad2
Input voltage	DC+12 V
Input current	110 mA
Max power consumption	1.32 W
Dimensions	W90 x D138 x H16 mm


For detailed installation instructions, refer to “4 Mounting Multi Pad2”.

3.2 NLZ BF

3.2.1 NLZ BF Control Board


Name	NLZ BF
Input voltage range	DC+32 V ±10 %
Input power	150 VA over
Max. output voltage	29 V peak
Rotation speed range	NLZ: 1,000-40,000 min ⁻¹ (General mode) 100-6,000 min ⁻¹ (Rotary Endo mode)
Torque range	NLZ: 4.2 Ncm (1:1) (General mode) 0.3-3.0 Ncm (1:1) (Rotary Endo mode)
Dimensions	W115 x D68.5 x H40 mm
External Electrical Power Source	IEC 60601-1 and Class II is required

For installing the NLZ BF Control Board with your dental chair table, use a screw or use Velcro (Magic Fastener).

3.2.2 AC – DC Adapter Board


Name	AC – DC Adapter Board
Input voltage range	AC 24 V ± 10 %
Input power	150 VA over
Max. output voltage	DC+32 V ± 10 %
Dimensions	W60.2 x D81.4 x H45.2 mm (Size of exclusive case)
External Electrical Power Source	IEC 60601–1 and Class II is required

For installing the AC – DC Adapter Board with your dental chair table, use the screws.

3.2.3 Selector board

Name	NLZ SEL
Input voltage range	DC+12 V
Input current	0.2 A
Dimensions	W68.5 x D56.5x H38.5 mm

4 Mounting Multi Pad2

4.1 Mounting MP2 STAY

⚠ CAUTION

- This product can be installed on the desired place. You need to drill holes on the installation place. Be sure to select a level, stable and safe surface to install the product.

- (1) Insert the post A into the hole of the mounting stay and tighten with the fixing nut.


Figure 1.

- (3) Align the screw hole of the clamp holder and post B and insert the screw with washer.


Figure 3.

- (2) Attach the stay ass'y with the screws and nuts at the bottom of the tray table or mobile cart unit.


Figure 2.

- (4) Tighten the screw with washer.
- (5) Attach the clamp holder to Multi Pad2 with two screws.


Figure 4.

⚠ CAUTION

- Be sure to use the two screws to fix the stay ass'y.

4 Mounting Multi Pad2

(6) Insert the post B into the stay ass'y.


Figure 5.

(8) Attach the cover onto the connection part of the post A and clamp holder.


Figure 7.

(7) Align the flat part of post B and the screw hole of post A and tighten the adjustment screw to fix the post B.


Figure 6.

(9) Bundle the extra cable with the attached cable tie, to complete the installation.


Figure 8.

⚠ CAUTION

- Check for loose nuts and screws and tighten them periodically.
- After installation is completed, check that the Stay and the Control Unit are stable and secure. If any abnormality is found, stop using the product and contact your Authorized NSK Dealer.

4.2 Mounting Angle Plate

- (1) Place the Angle Plate aligning with the screw holes of Multi Pad2.
- (2) Tighten four attached screws.
- (3) Mount the Angle Plate to the dental chair table.


Figure 9.

Dimensions of Angle Plate

(mm)


Figure 10.

4.3 Mounting Ferrite Core

- (1) Mount the Ferrite Core to the NLZ BF cables with one turn.
- (2) Place the UNITY BAND into the slot on the side of the Ferrite Core, and tie the cables.
- (3) Trim off and discard the excess length from the UNITY BAND.


Figure 11.

5 Connection

The NLZ BF controls the motor speed either by air pressure or external voltage.

5.1 Multi Pad2 and NLZ BF


Figure 12.

5.1.1 Connecting Harness Assy to NLZ BF

Connect the NLZ POWER HARNESS attached with the Multi Pad2 to the NLZ BF (Figure 13).


Figure 13.

5.1.2 DIP switch setting

NLZ BF has a 8 bits DIP switch.

Bit3 switch is for selection for motor speed control, either by external voltage or air pressure from the dental chair table.

For detailed settings for speed control, refer to “5.1.3 NLZ BF speed control settings”.

DIP switches

Specifications can be changed by changing the DIP switch settings.


Figure 14.

Bit	Description	Function	Details
1	Not used	Always keep turned OFF	Operation is not guaranteed if turned ON
2	Air drive	ON: Air drive OFF: Remote / potentiometer drive	Turn ON if driving motor according to air pressure
3	LED presence / absence	ON: LED absent OFF: LED present	Turn ON if not using LED in NLZ motor
4	Not used	Always keep turned OFF	Operation is not guaranteed if turned ON
5	Not used	Always keep turned OFF	Operation is not guaranteed if turned ON
6	Not used	Always keep turned OFF	Operation is not guaranteed if turned ON
7	Number of motors	ON: Two motors OFF: One motor	Turn ON if using two motors
8	Not used	Always keep turned OFF	Operation is not guaranteed if turned ON

5.1.3 NLZ BF speed control settings

The NLZ BF has two selection of motor speed control either by an air pressure or an external voltage.

Air pressure control

When you control NLZ's motor speed by air, it is necessary to set Bit2 of the DIP switch to ON. Then connect air tube from Drive air to air pressure sensor (IC302) on the NLZ BF as shown below.



Figure 15.

Put the Metal Tube Holder (Tube Retainer) through tube (left side in Figure 16). Make sure the holder is put in the direction shown below (right side in Figure 16).

Then, mount the tube onto the sensor. After mounting the tube, move the tube holder over the sensor for secure holding.

If the connection is loose, there is a risk of air leakage.

Figure 16.

5 Connection

If you need to use either Ø2.5 or Ø1.5 tube, use the tube relay (Mini Fitting) included in the product (Figure 17).

Figure 17.

Tube Holder, Mini Fitting and Ø1.5 tube are included in NLZ BF SET.

■ Air pressure sensor

The default setting for air pressure and motor rotation speed is shown as the figure below (Figure 18). This setting can be changed with the Multi Pad2.

For the air pressure adjustment, refer to the Multi Pad2 operation manual “6.1 Setting the supply pressure value with pneumatic dental chair unit with foot control”.


Figure 18.

⚠ CAUTION

- This product is an air pressure controlled electric motor system. Adjust the air pressure of the dental chair table to 0.392 MPa (4 kg/cm²) or less. Always check moisture condensation from the dental chair table. If too much moisture in the air, the system may not deliver the designed performance. Remove the moisture.

External voltage control

NLZ's motor can be controlled without using Multi Pad2 by connecting the toggle switches and potentiometer to EX-IF HARNESS (U1104501). In this case, only one motor of NLZ 1 in General mode can be used. Refer to Figure 31 in “9 Wiring” for the wiring method.


Figure 19.

5.2 Multi Pad2


Figure 20.

When only one switch is closed, that tool is selected.
If two or more tools are closed at the same time, a warning will be displayed on the LCD of Multi Pad2.

Color	Signal name	Function	Remarks
GREEN	+5 V	DC 5 V	Do not use this cable for the purpose other than potentiometer connection. Not used for air pressure control.
PINK	Speed_VR	Motor speed adjusting voltage input terminal	Variable within arrange of DC 0 V to +5 V (Potentiometer: 5 kΩ (B curve is recommended)) Not used for air pressure control.

5 Connection

Color	Signal name	Function	Remarks
WHITE	NLZ1	Motor select (NLZ1)	SW or Transistor connect
BLACK	NLZ2	Motor select (NLZ2)	SW or Transistor connect
GRAY	-	-	Not used
BLUE	GND	GND	-


Figure 21.

6 Power Supply Connection

The power of Multi Pad2 is sourced from the NLZ BF.

Power supply connection to the NLZ BF

Use CN100 for connection with the AC-DC Adapter Board CN401.

When you use DC32 V for the power source, do not use the AC-DC Adapter Board, and connect the power source to the connector (CN100) directly. We recommend to use U1104504 Fuse Harness for connecting CN401 on the AC-DC adapter board and NLZ BF. U1104504 Fuse Harness has a 6.3 A fuse (Littelfuse 021506.3MXP).

Precaution Notice

- Cord connection must be done while the power is off.
- A surge may flow into the fuse which may result in the damages on the PC boards if the connection is made while the power is on.
- Select the right insulated transformer and the switching regulator in accordance with the IEC 60601-1 standard because of no insulation of primary side in the AC-DC circuit.

CN400 (AC_IN)

No.	Signal name	Function	Remarks
1	AC24_1	AC 24 V input terminal	AC 24 V $\pm$ 10 % / 150 VA over
2	-	Not used	-
3	AC24_2	AC 24 V input terminal	-

AC-DC Adapter Board CN401 (DC OUT)

No.	Signal name	Function	Remarks
1	DC32 V	DC32 V output terminal	-
2	GND	Ground terminal	-

7 Cautions When Installing to Delivery Table

NLZ BF Control Board

- Ensure correct location for installation of the NLZ BF. Avoid water and exposure to heat.
- Ensure mounting space is dry and not exposed to water or oil.
- Always keep the label side of the NLZ BF to the top, or 90 degrees angle for air ventilation. Do not keep the label side of the NLZ BF to the bottom or more than 90 degrees angle.
- The power supply must be met the IEC 60601-1 requirement. Use the power supply for the double insulation type.

AC-DC Adapter Board

- Make sure to select the right spot to install the AC-DC Adapter Board. Avoid the AC-DC Adapter Board for water and space exposing heat.

Multi Pad2

- Ensure correct location for installation of the Multi Pad2. Do not install the Multi Pad2 near the chair / unit handle.
- Make a necessary angle adjustment of the Multi Pad2 for improved visibility.
- When connecting cord, allow extra length. Do not pull the cord tightly.

Motor cord

- When connecting cord, allow extra length (service loop). Do not pull the cord tightly.
- Fix the cord securely by using cord holder / fixer inside the delivery table.

8 EMC Information (Electromagnetic Compatibility Information)

The environments of use: Professional healthcare facility environment

Guidance and manufacturer's declaration - Electromagnetic Emissions.		
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment. Guidance and manufacturer's declaration electromagnetic emissions.		
Emissions test	Compliance	Electromagnetic environment - guidance
Radiated emissions CISPR 11/EN 55011	Group 1 Class B	The product uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emissions CISPR 11/EN 55011	Group 1 Class B	The product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2/EN 61000-3-2	Class A (220-240 V)	
Voltage fluctuations/flicker emissions IEC 61000-3-3/EN 61000-3-3	Complies (220-240 V)	

Guidance and manufacturer's declaration - Electromagnetic Immunity		
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.		
Immunity test	IEC/EN 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2/EN 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4/EN 61000-4-4	Input a.c. power PORT ± 2 kV Signal input/output parts PORT ± 1 kV	Input a.c. power PORT ± 2 kV Signal input/output parts PORT ± 1 kV
Surge IEC 61000-4-5/EN 61000-4-5	Line to line ± 0.5 kV, ± 1 kV Line to ground ± 0.5 kV, ± 1 kV, ± 2 kV	Line to line ± 0.5 kV, ± 1 kV Line to ground ± 0.5 kV, ± 1 kV, ± 2 kV

8 EMC Information (Electromagnetic Compatibility Information)

Voltage dips IEC 61000-4-11/EN 61000-4-11	0 % U_T ; 0.5 cycle (0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°) 0 % U_T ; 1 cycle and 70 % U_T ; 25 (50 Hz) /30 (60 Hz) cycles Single phase ; 0°	0 % U_T ; 0.5 cycle (0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°) 0 % U_T ; 1 cycle and 70 % U_T ; 25 (50 Hz) /30 (60 Hz) cycles Single phase ; 0°
Voltage interruptions IEC 61000-4-11/EN 61000-4-11	0 % U_T ; 250 (50 Hz) /300 (60 Hz) cycle	0 % U_T ; 250 (50 Hz) /300 (60 Hz) cycle
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8/EN 61000-4-8	30 A/m	30 A/m
Proximity magnetic field IEC 61000-4-39	134.2 kHz (Pulse modulation) 2.1 kHz, 50 % Duty, 13.56 MHz (Pulse modulation) 50 kHz 50 % Duty	134.2 kHz (Pulse modulation) 2.1 kHz, 50 % Duty, 13.56 MHz (Pulse modulation) 50 kHz 50 % Duty
NOTE U_t is the a.c. mains voltage prior to application of the test level.		

8 EMC Information (Electromagnetic Compatibility Information)

Guidance and manufacturer's declaration - Electromagnetic Immunity			
The product is intended for use in the electromagnetic environment specified below.			
The customer or the user of the product should assure that it is used in such an environment.			
Immunity test	IEC/EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6/ EN 61000-4-6	3 V rms 0.15 MHz - 80 MHz 6 V rms ISM bands between 0.15 MHz and 80 MHz	3 V rms 6 V rms	WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the product, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3/ EN 61000-4-3	3 V/m 80 MHz - 2.7 GHz 9 V/m in ISM bands 710 MHz, 745 MHz, 780 MHz, 5.24 GHz, 5.5 GHz, 5.785 GHz 27 V/m in ISM band 385 MHz 28 V/m in ISM bands 450 MHz, 810 MHz, 870 MHz, 930 MHz, 1.72 GHz, 1.845 GHz, 1.97 GHz, 2.45 GHz	3 V/m 9 V/m 27 V/m 28 V/m	
NOTE 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the product series is used exceeds the applicable RF compliance level stated above, the product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product series.		
b	Over the 150 kHz to 80 MHz frequency range, the field strength should be less than 3 V/m.		

Cables and accessories	Maximum length	Complies with
Micromotor with Motor Cord	2.2 m	RF emissions, CISPR 11: Class B/ Group 1

9 Wiring

9.1 Wiring Diagram

NLZ Built-in Motor System Diagram (2motors) Multi Pad2 Control



Figure 22.

*1 Refer to "9.2 Connecting joint".

Unnecessary wires should be pulled out from the connector or insulated and bundled so as not to short-circuit.

Figure 23.

NLZ Built-in Motor System Diagram (2motors) Air Control




Figure 24.

*1 Refer to "9.2 Connecting joint".

Unnecessary wires should be pulled out from the connector or insulated and bundled so as not to short-circuit.

MODEL
(ORDER CODE)

*1 Refer to "9.2 Connecting joint".

Unnecessary wires should be pulled out from the connector or insulated and bundled so as not to short-circuit.

NLZ Built-in Motor System (1motor) Air Control

AC-DC adapter board
NLZ ADP (U1104051) w/ Case

Not required for dental chairs
that can supply DC+32V±10%.

FUSE HARNESS
(U1104504)

Power Cable w/ Fuse

Fuse : 021506.3MXP
(Littelfuse)

Double check the
right polarity.


Connect the air tube from drive air to IC302 (Air
pressure sensor)

- DIP switch: Bit2 ON

- Drive air pressure: Max. 0.4 Mpa / 4.0 bar / 58.02 psi

- Nozzle diameter of Air pressure sensor: 3mm

AIR INPUT

OVER AWG20

Input AC24V±10%
(Over 150VA)


NLZ POWER
HARNESS
(U1104500)

Black (GND)
Red (DC+32V)


GND DC+32V

Control Board
NLZ BF (U1104002) w/ Case

Dip switch

Multi Pad2
(UA1408-0001)


Ferrite Core
(D4Z400201)


Not used
Purple (GND)
Yellow (RIN)
Orange (TOUT)

Tool Holder for Auto Selection

Open : Not selected
Close : Selected


MODEL
(ORDER CODE)

*1 Refer to "9.2 Connecting joint".

Unnecessary wires should be pulled out from the connector or insulated and bundled so as not to short-circuit.

MODEL
(ORDER CODE)

When only one switch is Close, that tool is selected.
If two or more tools are closed at the same time, a warning will be displayed on the LCD of Multi Pad2.

Unnecessary wires should be insulated and bundled so as not to short-circuit.


Figure 30.

NLZ Built-in Motor System Diagram (1motor) External Interface Control

MODEL
(ORDER CODE)

AC-DC adapter board
NLZ ADP (U1104051) w/ Case

Not required for dental chairs
that can supply DC+32V±10%.

*1 Refer to "9.2 Connecting joint".

Unnecessary wires should be pulled out from the connector or insulated and bundled so as not to short-circuit.

9.2 Connecting joint

If the motor cord tube and dental chair tube have different diameters, mount the joint as shown in the figure below.


Figure 32.

⚠ CAUTION

- When connecting the tube to the joint, insert it firmly to the end. Incomplete insertion may cause water or air leakage.
- Do not repeatedly insert and remove the tube. Doing so may cause loosening or damage to the joints, resulting in water or air leakage.

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