OM-DE0907EN 000

Micromotor with LED

NLZ

OPERATION MANUAL

Rx Only

MADE IN JAPAN

1 Preface

1-1 User and indications for use

■ User

Dentist, dental hygienist

■ Indications for Use

The NLZ Built-in Motor System is indicated for use in the field of prophylaxis dentistry, restorative applications including cavity preparation and endodontic therapy, prosthodontics applications such as crown preparations.

1-2 Related manuals

Operation Manual Title	Operation Manual No.
Multi-Control Panel for built-in system Multi Pad2 OPERATION MANUAL	OM-DE0989EN
Multi Pad2 NLZ Built-in Motor System INSTALLATION GUIDE	OM-TM0903EN

See the “1-4 Symbol” section on page 5-6 for reference.

1–3 Precautions for handling and operation

- Read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
⚠ WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.

⚠ CAUTION

Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.

NOTICE

General product specification information highlighted to avoid product malfunction and performance reduction.

⚠ WARNING

- The motor is delivered in a non-sterile condition and must be cleaned and sterilized prior to its first use and after each patient use.
- Do not operate this product close to patients with cardiac pacemakers as there is a danger that it may affect the pacemaker.

⚠ CAUTION

- When operating this motor always consider the safety of the patient.
- The motor is designed only for clinical dental use by qualified personnel.
- This motor cannot operate without exclusive control unit.
- Do not allow any impact to the product. Do not drop the product. Doing so could cause deformation or malfunctions.
- Make sure that each part is properly connected. Improper connections might cause faulty operation and water or air leakage.

⚠ CAUTION

- Prior to clinical use, inspect the motor. Check for vibration, noise and overheating. If any abnormalities are found, stop using the motor immediately and contact your authorized NSK dealer. (Refer to “4 Check before Operation”)
- Do not allow water to enter the motor. It may cause malfunction of the motor.
- Do not exceed the drive motor speed recommended by the handpiece (option) manufacturers.
- Do not use air contaminated by dust, moisture or oil. If the air contains water or dust, it might cause a malfunction or overheating.
- Do not use this motor under strong stress for long time. It causes overheating.

CAUTION

- Do not sterilize the motor cord.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Surface temperature shall be more than 60 °C (140 °F) without using coolant air, to avoid this event, be sure to use coolant air.
- Follow the instruction of the dental chair unit manufacturer for cleaning and disinfection of the water line inside the motor cord (DUWL).
- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information (refer to “10 EMC Information (Electromagnetic Compatibility Information)”).

CAUTION

- This product may be affected when used in rooms with electromagnetic interference and portable or mobile RF communication equipment. Do not use this product around any equipment that generates electromagnetic waves such as an ultrasonic generator and electrical scalpel. When using equipment generating electromagnetic waves, turn off the power to the dental chair unit (refer to “10 EMC Information (Electromagnetic Compatibility Information)”).

CAUTION

- The use of accessories other than those specified, with the exception of accessories sold by the manufacturer of this product as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of this product.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.
- Please report any serious incident that may occur in connection with this device and its accessories to the manufacturer as well as to your national authority.

1-4 Symbol

Wear protective gloves

Wear a face mask

Wear eye protection

This product can be sterilized in a steam sterilizer at 135 °C

Consult operation instructions

Caution

Marking on the outside of Equipment or Equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment

Date of manufacture

Conforms to CE
European Directive of Medical device
directive 93/42/EEC

Distributor

Manufacturer

Medical Device

Caution: U.S. Federal law restricts this device to be used by or on the order of a licensed dental professional

Device Name for Health Canada

Catalog number (Order Code)

Serial No.

Unique Device Identifier

GS1 DataMatrix for Unique Device Identifier

Temperature limitation

Humidity limitation

Atmospheric pressure limitation

1-5 Disposing of products

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

1-6 Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non-NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model. Contact your authorized NSK dealer if repairs are necessary.

2 Product Description

2-1 Package contents

No.	Part Name	Quantity	Remarks
1	O-ring	1	Blue
2	O-ring	3	Black
3	O-ring	5	Black

*Handpiece is not included.

2-2 Applicable product

Product Name	REF	Remarks
NLZ CDB	EA1152X061	Motor Hose
NLZ Endo	C1130	Handpiece

3 Preparation for Use

3-1 Connecting and disconnecting the motor and motor cord

- Align and insert firmly the motor pin into the pin holes of the motor cord connector, and fasten the motor cord nut securely. (Figure 1.)


Figure 1.

- 2 Make sure the motor is connected firmly to the motor cord.
- 3 Rotate the motor and check for vibration, noise and overheating.

To disconnect, unscrew and detach the motor cord nut, and gently pull out the motor cord.

⚠ CAUTION

- Do not connect or disconnect the motor until the motor has completely stopped.
- Do not depress the foot pedal when disconnecting the motor.

NOTICE

- The NLZ motor is intended for use with connecting to NLZ joint (E1152701).

3-2 Using barrier sleeve (motor cord)

Use a disposable barrier sleeve that meets the U.S. and Canada regulation for the motor and motor cord, in order to prevent cross contamination between each patient. Install the disposable barrier sleeve to the motor (including the motor cord) then fix the barrier sleeve with a tape. (Figure 2.)


Figure 2.

A barrier sleeve with a diameter of larger than the diameter of the motor ($\varnothing 22$ mm) is recommended.

⚠ CAUTION

- Dispose of the barrier sleeve after each patient.
- For details of the barrier sleeve, confirm the operation manual issued by manufacturer.

3-3 Connecting and disconnecting the handpiece from the motor

- 1 Insert the E-type handpiece into the motor insert, and turn until it clicks as the positioning pin on the handpiece falls into positioning hole on the motor. (Figure 3.)


Figure 3.

- 2 Confirm that the handpiece is securely connected to the motor.

- 3 Rotate the motor and check for vibration, noise and overheating.

To disconnect, hold the motor and the handpiece, then pull apart.

CAUTION

- Care should be taken when using any handpiece attachment which requires lubrication. Following lubrication stand the handpiece attachment vertically until the excess solution drains. Wipe the handpiece then attach it to the motor.
- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Do not rotate the motor when connecting or disconnecting the handpiece.

4 Check before Operation

Follow the check below before use. If any abnormalities are found, stop using the motor immediately and contact your authorized NSK dealer.

- 1 Connect the handpiece and make sure that it is securely connected.
- 2 Rotate the motor for about one minute at the Max. rotation speed. During rotation, check for abnormalities such as abnormal rotation, vibration and noise.
- 3 After the motor rotation has completely stopped, touch the motor to confirm there is not heating abnormally.

5 Post-use Maintenance

After each patient, maintain the product as follows.

WARNING

- To ensure cleaning and sterilization efficacy, use only the following procedures for reprocessing.
- For cleaning and disinfection, use disinfectant ethanol or disinfectant isopropanol that does not contain any additives (hereinafter referred to as disinfectant alcohol). Use of disinfectants other than disinfectant alcohol may cause discoloration or cracking.
- Follow any additional local directives, standards, and guidelines for cleaning and sterilization.

WARNING

- Clean the motor immediately (within 1 hour) after each treatment to remove residue. Failure to properly maintain the motor may cause infection, product failure, overheating leading to burn injuries.

CAUTION

- Do not use the following liquids to wipe, immerse, or clean the product: strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.
- Do not immerse NSK instruments in disinfectant solution or clean in ultrasonic devices.

⚠ CAUTION

- Ensure that no liquid enters the motor during cleaning and disinfecting.
- Do not use a washer-disinfector to clean the motor.
- Do not lubricate the motor.
- Use appropriate disinfectant alcohol according to the regulations and guidelines of each country.

5-1 Preparation for maintenance EN

Items to prepare

- Eye protection
- Mask
- Gloves
- Cloth dampened with disinfectant alcohol
- FDA-approved sterilization pouch that conforms to ISO 11607-1

- 1** Wear eye protection, a mask and gloves to prevent infection (Figure 4.).

Figure 4.

- 2 Remove the bur (Figure 5.).


Figure 5.

- 3 Wipe the exterior of the handpiece and motor clean using a cloth dampened with disinfectant alcohol (Figure 6.). When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.


Figure 6.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties. The following disinfectants can be used in the United States and Canada.

Example) CaviCide[®], CaviWipes[®] (manufactured by Metrex)

CaviCide[®] and CaviWipes[®] are registered trademarks of Metrex Research, LLC.

- 4 Operate the motor at the chair-side for at least 20 seconds to purge fluids in the motor (Figure 7.).

Figure 7.

- 5 Remove the handpiece from the motor.
For reprocessing of the handpiece, refer to the operation manual of the handpiece's manufacturer.
- 6 Remove the barrier sleeve.
- 7 Remove the motor from the motor cord (Figure 8.).

Figure 8.

⚠ CAUTION

- Do not twist, kink, or squeeze the motor cord. Doing so may cause malfunction.

5-2 Cleaning

Wipe the exterior of the motor clean using a cloth dampened with disinfectant alcohol (Figure 9.). When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.


Figure 9.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties. The following disinfectants can be used in the United States and Canada.

Example) CaviCide[®], CaviWipes[®] (manufactured by Metrex)
CaviCide[®] and CaviWipes[®] are registered trademarks of Metrex Research, LLC.

5-3 Sterilization

- 1 Insert the motor into an FDA-approved sterilization pouch that conforms to ISO 11607-1, and seal the pouch.
- 2 Perform steam sterilization with the following conditions.

Type	Gravity Displacement	Pre-Vacuum (Dynamic Air Removal)
Temperature	132 °C (270 °F)	132 °C (270 °F)
Full Cycle Time	15 min	4 min
Drying Time	30 min	30 min

- 3 Allow the load to dry for approximately 30 minutes or more before removing from the sterilizer. Ensure the load is completely dry and no water droplets are visible on the sterilization pouch prior to removal from the sterilizer.

* For more information about sterilization standards, refer to the medical device regulations specific to your location.

CAUTION

- Use an FDA-approved steam sterilizer to perform sterilization.
- Follow local rules, regulations, and guidelines regarding the reprocessing of devices.
- Do not use sterilization pouches containing water-soluble adhesive components such as PVA (polyvinyl alcohol). Adhesive components washed out during sterilization may enter the product, causing problems such as poor rotation, sticking and failure to operate correctly.

⚠ CAUTION

- Do not touch the product immediately after sterilization as it will be very hot and must remain in a sterile condition.
- Do not perform steam sterilization on the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Clean and disinfect before you sterilize the motor to ensure the sterilization process is effective and to prevent damage to the motor.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- Do not sterilize the motor cord.

⚠ CAUTION

- Be sure to use sterilizers that can perform sterilization up to 135 °C (275 °F). In some sterilizers, the chamber temperature may exceed 135 °C (275 °F). Do not use these sterilizers as failure of the motor could occur. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Steam sterilization is recommended for the product. The validity of other sterilization methods (such as plasma sterilization or EOG sterilization) is not confirmed.
- Handle the motor with care and avoid scratching the optic outlet.

NOTICE

- NSK recommends sterilization according to EN 13060 (class B). Other sterilization methods may reduce the life span of your motor.
- Follow the manufacturer's Instructions for Use that came with your sterilizer.

5-4 Storage

Store the product in a dry, clean location.

⚠ CAUTION

- After the sterilization and drying cycles are complete, remove the motor immediately from the sterilizer to store it.

⚠ CAUTION

- Store the product in a well ventilated place out of direct sunlight and within the range of temperature, humidity and pressure specified in "9 Product Specifications". The air should be free from dust, salt and sulphur.
- Sterilization is not guaranteed after the sterilization retention period specified by the manufacturer and seller of the sterilization pouch has elapsed. If the sterilization retention period has elapsed, perform sterilization again with a new sterilization pouch.

6 Maintenance

6-1 Replace the O-rings (motor insert)

If a deteriorated O-ring makes it difficult to attach the handpiece or causes water or air to leak, replace the O-ring. Using a needle or other pointed tool, remove defective O-rings from the motor insert section, and fit new O-rings into the O-ring grooves (Figure 10.).


Figure 10.

⚠ CAUTION

- Do not force the new replacement O-ring with excessive pressure.
- The blue O-ring is thinner than the three black O-rings. When inserting new O-rings, make sure they are inserted in the correct grooves as shown in Figure 10.

6-2 Replace the O-rings (motor cord connection)

If you detect an air/water leak coming from the motor cord connection, replace the O-rings.

Remove the damaged O-rings in the rear of motor with a pointed tool and insert new O-rings into the O-ring grooves. (Air circuit: 2 pcs., water circuit: 2 pcs., coolant air circuit: 1 pc.) (Figure 11.)


Figure 11.

6-3 Periodical maintenance checks

Perform periodical maintenance checks every three months, referring to the check sheet below. If any abnormalities are found, contact your authorized NSK dealer.

Points to check	Details
Connection	Check the connection between the motor and the handpiece for any problems.
Rotation	Rotate the motor and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.

7 Troubleshooting

When a problem is detected, check the following again before requesting a repair. If none of these is applicable or if the trouble is not remedied even after an action has been taken, a failure of this product is suspected. Contact your authorized NSK dealer.

Problem	Cause of problem	Action
Motor does not run.	Motor cord is not connected correctly.	Check the connection.
LED does not light.	Reached the end of life expectancy.	Contact your dealer.
Air/Water leakage.	Motor cord is not connected correctly.	Check the connection.
	O-ring has been damaged.	Replace with new O-ring.
The motor heats abnormally.	Air pressure for cooling from the dental chair unit is not applied or insufficient.	Check air pressure of the dental chair unit.

* Repair cannot be made by the customer.

8 Spare Parts List

Model	REF	Remarks
O-ring Set	E1135083	For motor insert Black: 3 pcs. / Blue: 1 pc.
O-ring	D0313084070	For motor insert Blue: 1 pc.
O-ring	D0312074080	For motor insert Black: 1 pc.
O-ring	D0312010050	For motor cord connection

9 Product Specifications

■ Specification table

Model	NLZ
Max. Rated input Electric voltage	AC 27 V
Max. Rated input Electric current	4.5 A
Max. Rotation Speed	40,000 min ⁻¹
Max. Torque	4.2 N·cm
Dimensions	Ø22xH68.8 mm
Weight	57 g
Spray Supplying Water	More than 65 mL/min (0.2 MPa)
Chip Air	More than 1.5 NL/min (0.2 MPa)
Coolant Air	More than 6.5 NL/min (0.25 MPa)
Operation Time	ON: 3 min, OFF: 10 min
Motor Connection Type	ISO3964: Type 3 (short)

■ Use, transportation and storage environment

Use Environment	Temperature: 10-40 °C (50-104 °F), Humidity: 30-75 % (No Condensation), Pressure: 700-1,060 hPa
Transportation and Store Environment	Temperature: -10-50 °C (14-122 °F), Humidity: 10-85 %, Pressure: 500-1,060 hPa

10 EMC Information (Electromagnetic Compatibility Information)

Guidance and manufacturer's declaration - Electromagnetic Emissions.		
The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that is used in such an environment. Guidance and manufacturer's declaration electromagnetic emissions.		
Emissions test	Compliance	Electromagnetic environment - guidance
Radiated emissions CISPR 11/EN 55011	Group 1 Class B	The product uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Conducted emissions CISPR 11/EN 55011	Group 1 Class B	
Harmonic emissions IEC 61000-3-2/EN 61000-3-2	Class A (220-240 V)	
Voltage fluctuations/flicker emissions IEC 61000-3-3/EN 61000-3-3	Complies (220-240 V)	

Guidance and manufacturer's declaration - Electromagnetic Immunity

The product is intended for use in the electromagnetic environment specified below.
The customer or the user of the product should assure that it is used in such an environment.

Immunity test	IEC/EN 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2/EN 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4/EN 61000-4-4	Input a.c. power PORT ± 2 kV Signal input/output parts PORT ± 1 kV	Input a.c. power PORT ± 2 kV Signal input/output parts PORT ± 1 kV
Surge IEC 61000-4-5/EN 61000-4-5	Line to line ± 0.5 kV, ± 1 kV Line to ground ± 0.5 kV, ± 1 kV, ± 2 kV	Line to line ± 0.5 kV, ± 1 kV Line to ground ± 0.5 kV, ± 1 kV, ± 2 kV

Voltage dips IEC 61000-4-11/EN 61000-4-11	0 % U_T ; 0.5 cycle (0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°) 0 % U_T ; 1 cycle and 70 % U_T ; 25 (50 Hz) /30 (60 Hz) cycles Single phase ; 0°	0 % U_T ; 0.5 cycle (0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°) 0 % U_T ; 1 cycle and 70 % U_T ; 25 (50 Hz) /30 (60 Hz) cycles Single phase ; 0°
Voltage interruptions IEC 61000-4-11/EN 61000-4-11	0 % U_T ; 250 (50 Hz) /300 (60 Hz) cycle	0 % U_T ; 250 (50 Hz) /300 (60 Hz) cycle
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8/EN 61000-4-8	30 A/m	30 A/m
Proximity magnetic field IEC 61000-4-39	134.2 kHz (Pulse modulation) 2.1 kHz, 50 % Duty, 13.56 MHz (Pulse modulation) 50 kHz 50 % Duty	134.2 kHz (Pulse modulation) 2.1 kHz, 50 % Duty, 13.56 MHz (Pulse modulation) 50 kHz 50 % Duty
NOTE U_t is the a.c. mains voltage prior to application of the test level.		

Guidance and manufacturer's declaration - Electromagnetic Immunity

The product is intended for use in the electromagnetic environment specified below.

The customer or the user of the product should assure that it is used in such an environment.

Immunity test	IEC/EN 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6/EN 61000-4-6	3 V rms 0.15 MHz - 80 MHz 6 V rms ISM bands between 0.15 MHz and 80 MHz	3 V rms 6 V rms	<p>WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the product, including cables specified by the manufacturer.</p> <p>Otherwise, degradation of the performance of this equipment could result.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3/EN 61000-4-3	3 V/m 80 MHz - 2.7 GHz 9 V/m in ISM bands 710 MHz, 745 MHz, 780 MHz, 5.24 GHz, 5.5 GHz, 5.785 GHz 27 V/m in ISM band 385 MHz	3 V/m 9 V/m 27 V/m	

	28 V/m in ISM bands 450 MHz, 810 MHz, 870 MHz, 930 MHz, 1.72 GHz, 1.845 GHz, 1.97 GHz, 2.45 GHz	28 V/m	
--	---	--------	--

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the product series is used exceeds the applicable RF compliance level stated above, the product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product series.

b Over the 150 kHz to 80 MHz frequency range, the field strength should be less than 3 V/m.

Cables and accessories	Maximum length	Complies with
Micromotor with Motor Cord	2.2 m	RF emissions, CISPR 11: Class B/ Group 1

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NSK America Corp. www.nskdental.com