iProphy Air

iProphy Air

Rx Only

Please read this Operation Manual carefully before use and file for future reference

MADE IN JAPAN (GM-M0129E) 001

1. User and Indications for Use

User: Qualified Professionals

Indications for Use: The Nakanishi iProphy Air is a prophylaxis handpiece for use with disposable prophy angles in hygiene operatory to perform cleaning and polishing procedures on teeth.

2. Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

WARNING	Degree of Risk
WARNING	A hazard that could result in serious injury or damage to the device if the safety instructions are not followed.
CAUTION	A hazard that could result in light or moderate injury or damage to the device if the safety instructions are not followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

WARNING

- Clean and lubricate the handpiece immediately (within 1 hour) after each treatment to remove residue. Failure to properly maintain the handpiece may cause infection, product failure, overheating leading to burn injuries. (Refer to "8. Post-use Maintenance")

CAUTION

- Read this Operation Manual before use to fully understand the product functions.
- When operating the handpiece always consider the safety of the patient.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Do not attempt to disassemble the handpiece or tamper with the mechanism except as recommend by NSK in this Operation Manual.
- Do not allow any impact on to the handpiece. Do not drop the handpiece.
- Should the handpiece function abnormally cease operation immediately and contact your Authorized NSK Dealer.
- Do not use high acid water or sterilizing solutions to wipe, immerse or clean the handpiece.
- The products are delivered in a non-sterile condition and must be autoclaved prior to use.
- Perform function and maintenance checks periodically.
- If the product is not used for a long period check it is functioning correctly before using on a patient.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during surgery.
- US Federal law restricts the sale of this device to or by order of a licensed physician.

3. Accessory List

No.	Part Name	Quantity
1	E-Type Spray Nozzle	1

4. Setting the Air Supply Pressure

Measure the supply pressure at the handpiece/hose connection point and set the pressure to 0.20 - 0.25MPa as specified on the specification table. (Fig. 1)


Fig. 1

WARNING

- Do not exceed 0.25MPa.

CAUTION

- Do not use air contaminated by dust, moisture or oil.

5. Connection & Disconnection of the Handpiece

5-1 Connection

- 1) Insert the handpiece correctly into the hose connector and tighten the hose nut. (Fig. 2)
- 2) Make sure the handpiece is connected firmly to the Hose.


Fig. 2

5-2 Disconnection

- 1) Loosen the Hose Nut and and remove from the Hose. (Fig. 2)

6. Mounting and Removal of the Disposable Prophy Angle

CAUTION

- Do not attempt to mount or remove the Disposable Prophy Angle until the motor has completely stopped.

NOTICE

- The disposable Prophy Angle's positioning mark depends on the manufacturer. It may be different from the position of Fig. 3.

6-1 Mounting

To mount the handpiece and the Disposable Prophy Angle, align the positioning marks of the handpiece and the Disposable Prophy Angle and insert firmly.

CAUTION

After the Disposable Prophy Angle is mounted check that the positioning mark of Disposable Prophy Angle and handpiece are properly aligned.


Fig. 3

6-2 Removal

To disconnect hold the Disposable Prophy Angle and the handpiece and pull out straight.

7. Check before Treatment

Follow the check below before use. If any abnormalities are found stop using the handpiece immediately and contact your NSK dealer.

- 1) Connect the Disposable Prophy Angle. (Do not exceed maximum rotation speed for any item.)
- 2) Check that there is no abnormal rotation, vibration, noise or overheating.
- 3) If there are no immediate apparent abnormalities continue to rotate the handpiece for about 1 minute. Then check if there is no abnormal heating on the outer surface.
- 4) Operate it if there are no abnormalities above.

8. Post-use Maintenance

8-1 Preparation

- 1) Remove the Disposable Prophy Angle from the handpiece.
- 2) Wipe the exterior of the handpiece clean using a cloth moistened with disinfectant.
- 3) Remove the handpiece tip (IPA-HP) from the motor (IPA). (Fig. 4)
- 4) Remove the motor (IPA) from the hose.


Fig. 4

8-2 Cleaning

- 1) Clean the external surfaces of the handpiece (IPA-HP / IPA) with a soft bristled brush under running water.
- 2) Use a dry cloth to completely remove any exterior moisture.
- 3) Wipe the exterior of the handpiece (IPA-HP / IPA) clean using a cloth moistened with disinfectant.

CAUTION

- To clean the handpiece never use any solvent such as benzene or thinner.

8-3 Lubrication

Apply NSK PANA SPRAY Plus after every use and/or before autoclaving.

- 1) Mount the tip nozzle into the spray can port. (Fig. 5)
- 2) Insert the tip nozzle into the Drive Air Port of the motor (IPA). Hold the motor and spray for approximately 2-3 seconds. Apply lubricant until it comes out from the exhaust air port (see Fig. 6) for at least 2 seconds.
- 3) Attach the Spray Nozzle into the Spray Port nozzle on the can. (Fig. 7)
- 4) Insert the Spray Nozzle into the rear of the hanpeace tip (IPA-HP). Hold the handpiece tip and spray 2-3 times for approximately 2-3 seconds each time. Apply lubricant until it expels from the handpiece tip head for at least 2 seconds. (Fig. 7)


Fig. 5


Fig. 6


Fig. 7

CAUTION

- Firmly hold the handpiece to prevent it slipping when spray pressure is applied.
- Hold the spray can upright.

8-4 Sterilization

Sterilize the product by autoclave sterilization.

- 1) Insert the handpiece into an FDA-approved sterilization pouch that conforms to ISO 11607-1, and seal the pouch.
- 2) Perform steam sterilization with the following conditions.

Type	Gravity Displacement	Pre-Vacuum (Dynamic Air Removal)
Temperature	132°C (270°F)	132°C (270°F)
Exposure Time	15min	4min
Drying Time	30min	30min

CAUTION

- Use an FDA-approved steam sterilizer to perform sterilization.
- Follow local rules, regulations, and guidelines regarding the reprocessing of devices.
- Do not touch the product immediately after steam sterilization as it will be very hot and must remain in a sterile condition.
- Do not perform steam sterilization on the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Clean and lubricate the handpiece prior to sterilization. If blood remains on the internal surface it can become clotted and cause product failure.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- Be sure to use sterilizers that can perform sterilization up to 135°C. In some sterilizers, the chamber temperature may exceed 135°C. Do not use these sterilizers as failure of the handpiece could occur. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Steam sterilization is recommended for the product. The validity of other sterilization methods (such as plasma sterilization or EOG sterilization) is not confirmed.

8-5 Storage

Store the product in a dry, clean location.

CAUTION

- After the sterilization and drying cycles are complete, remove the handpiece immediately from the sterilizer to store it.
- Store the product in a well ventilated place out of direct sunlight and within the range of temperature, humidity and pressure specified in "10. Specifications".
- Sterilization is not guaranteed after the sterilization retention period specified by the manufacturer and seller of the sterilization pouch has elapsed.
- If the sterilization retention period has elapsed, perform sterilization again with a new sterilization pouch.


9. Periodic Maintenance Checks

Perform periodic maintenance checks every three months by referring to the check sheet below. If any abnormalities are found contact your Authorized NSK Dealer.

Points to check	Details
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, or overheating.

10. Specifications

Model	iProphy Air	
Hose Connection Type	ISO9168 Type 2	
Max Rotation Speed	2,000 - 3,000min ⁻¹ (0.25MPa)	
Drive Air Pressure	0.20 - 0.25MPa (2.0 - 2.5kgf/cm ²)	
Max Air Consumption	<66NL/min (0.25MPa)	
Use Environment	Temperature	10 - 35 °C
	Humidity	30 - 75%
Transportation and Store Environment	Atmospheric Pressure	700 - 1,060hPa
	Temperature	-10 - 50 °C
Environment	Humidity	10 - 85%
	Atmospheric Pressure	500 - 1,060hPa

11. Symbols

	This product is Autoclavable up to Max.135°C.
	Manufacturer
	Caution: US Federal law restricts the sale of this device to or by order of a licensed physician.

12. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the handpiece be used incorrectly or for anything other than the intended purpose, has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years after discontinuation of the model.

13. Option Parts List

Model	Order Code
MG-4H Multi Gauge	Z109400

14. Spare Parts List

Model	Order Code
iPA	M1020006
iPA-HP	M1020005
iProphy Air Color Ring Set	Y1002594
E-Type Spray Nozzle	Z019090

15. Disposing product

To avoid health risks and environmental contamination from the disposal of medical equipment a licensed physician must verify equipment is sterile prior to disposal. Seek the assistance of specialist firms who are licensed to dispose of specially controlled industrial wastes to dispose the product for you.

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