The NSK logo is displayed in a bold, blue, italicized sans-serif font.The text "Rx Only" is enclosed in a rounded rectangular box with a light blue gradient background.

Please read this Operation Manual carefully before use, and file for future reference.

Surgical Handpiece

Ti-Max X	X-SG20L
	X-SG25L
	X-SG93L
	X-SG93
S-Max	SG20

MADE IN JAPAN

OM-C0560E

Powerful Partners®

 **CAUTION for X-SG20L / SG20**

- This handpiece is intended for use only with Surgical Units with Torque Calibration System such as NSK Surgic Pro. Using this product with other surgical units may damage the handpiece or bone.

1. User and Indications for Use

User: Qualified Professionals

Indications for Use: Surgical use

2. Precautions for handling and operation

■ Please read these precautions carefully and use only as intended or instructed.

■ Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

 **WARNING**

- Immediately after a surgery, the handpiece should be cleaned, lubricated and sterilized. Failure to properly maintain the handpiece can lead to contaminants remaining inside the handpiece. Unremoved blood may result in product failure due to the internal blood coagulation. Improper maintenance may also lead to overheating, causing burn injuries, product failure, or difficulty in attaching and detaching the bur/drill. (Refer to "8. Post-use Maintenance")
- Depressing the Push Button while the handpiece is in rotation may lead to overheating, causing burn injuries or product failure. Avoid the push button to contact with any oral tissue.
- Supply coolant water while using the handpiece. No supplying the coolant water may lead to overheating, causing burn injuries or product failure.
- Do not use the handpiece for polishing purposes. Polishing paste could enter the handpiece, causing the failure of the Push Button or burn injuries due to the overheating of the handpiece head.
- Keep any debris or other foreign materials away from inside the gear or the handpiece. Foreign materials remaining inside may lead to overheating, causing burn injuries or other accidents.

 **WARNING**

- Ball bearings are wear items and may require replacement. Be sure to inspect the handpiece before use (Refer to "7. Checking the Handpiece Before Each Use"). If abnormal vibration or noise are found, the bearings or other internal parts may be worn out or damaged. Damaged parts may lead to handpiece overheating causing burn injuries. If overheating is suspected, stop using the handpiece immediately and contact your Authorized NSK Dealer.

CAUTION

- When operating the product always consider the safety of the patient.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Prior to clinical use, inspect the handpiece. Check for vibration, noise and overheating. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer. (Refer to "7. Checking the Handpiece Before Each Use")
- Do not disassemble or alter the handpiece except as recommended by NSK in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product.
- Operators and all others in the area must wear eye protection and a mask when operating this handpiece.
- Only X-SG20L / SG20 can use drill.
- Do not use bur/drills with problems listed below as the bur/drill may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur/drill.
 - Bur/drill which is cracked on the edge or axis.
 - Non-ISO(EN ISO) standard, or tampered bur/drill.
- Do not use bur/drills other than specified in "10. Specifications" as such bur/drills may be accidentally released or break during rotation.

CAUTION

- Always keep the bur/drill shank clean. Dirt or debris in the chuck could cause poor bur/drill concentricity or low chuck retention force.
- Always insert the bur/drill all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur may occur.
- Do not exceed Max. Bur/Drill Length shown in "10. Specifications".
- Always follow the instructions provided by the bur/drill manufacturer.
- Do not exceed the bur/drill speed recommended by the bur/drill manufacturer.
- Should the handpiece function abnormally during use, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.
- The handpiece is delivered in a non-sterile condition and must be autoclaved prior to use.
- Perform regular function and maintenance checks. (Refer to "9. Periodical Maintenance Checks")
- If the handpiece has not been used for a long period, rotate the handpiece and check for noise, vibration and overheating before use.

CAUTION

- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.
- The operation of the handpiece is permitted only on dental units which correspond to the standards IEC 60601-1 (EN 60601-1) and IEC 60601-1-2 (EN 60601-1-2).
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

X-SG93L / X-SG93

- Do not use a surgical bur with its maximum length of 25mm or longer. Using such bur with this handpiece may lead to premature abrasion of the bearings. In addition, sudden release, bend or breakage of the bur may occur.

3. Accessory List

No.	Part Name	Quantity	Remarks
1	Internal Spray Nozzle	1	Included in the package for X-SG20L / SG20
2	Nozzle Holder	1	Included in the package for X-SG20L / SG20
3	Y-Connector	1	Included in the package for X-SG20L / SG20
4	Cleaning Wire	1	Included in the package for X-SG20L / X-SG25L / SG20
5	E Type Spray Nozzle	1	-

4. Connecting and Disconnecting the Handpiece and the Motor

4-1 Connecting

- 1) Insert the handpiece direct to the motor (Fig. 1). Screw the handpiece until it clicks and locks into position.
- 2) Confirm that the handpiece is firmly connected to the motor.

4-2 Disconnecting

Hold the motor and the handpiece, then pull apart.


Fig. 1

CAUTION

- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Connect ONLY to E type motors (ISO 3964 (EN ISO 3964)).
- Do not exceed the Max. Rotation Speed (Motor) shown in "10. Specifications".

 **CAUTION**

X-SG20L / SG20

- This handpiece is intended for use only with Surgical Units with Torque Calibration System such as NSK Surgic Pro. Using this product with other surgical units may damage the handpiece or bone.

SG20

- The handpiece is designed for connection to a NSK Non-optic motor or to a motor with a stopper pin. Connection to any other type of motor cannot lock the handpiece.

5. Mounting and Removing the Bur/Drill

5-1 Mounting the Bur/Drill

X-SG20L / X-SG25L / SG20 (Fig. 2)

- 1) Insert the bur/drill until it is correctly seated in place.
- 2) Depress the Push Button and insert the bur into the chuck until the bur "notch" mechanism engages. Release the button.
- 3) Confirm that the bur/drill is not removed by gently pulling and pushing the bur/drill without depressing the Push Button.


Fig. 2

⚠ CAUTION

- Only X-SG20L / SG20 can use drill.

X-SG93L / X-SG93 (Fig. 3)

- 1) Insert the bur into the chuck.
- 2) Depress the Push Button to open the chuck (❶).
- 3) Insert the bur fully into the chuck until it stops (❷) then release the Push Button.
- 4) Confirm that the bur is not removed by gently pulling and pushing the bur without depressing the Push Button.


Fig. 3

NOTICE

- Grip the handpiece while placing your thumb tip on the Push button which makes it easier to depress the button.

5-2 Removing the Bur

Depress the Push Button to open the chuck (❶) and remove the bur/drill (❷).


Fig. 4

CAUTION

- Do not use bur/drills with problems listed below as the bur/drill may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur/drill.
 - Bur/drill which is cracked on the edge or axis.
 - Non-ISO(EN ISO) standard, or tampered bur/drill.
- Do not use bur/drills other than specified in "10. Specifications" as such bur/drills may be accidentally released or break during rotation.
- Always keep the bur/drill shank clean. Dirt or debris in the chuck could cause poor bur/drill concentricity or low chuck retention force.
- Always insert the bur/drill all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur may occur.
- Do not exceed Max. Bur/Drill Length shown in "10. Specifications".
- Do not exceed the bur/drill speed recommended by the bur/drill manufacturer.
- Do not mount or remove the bur/drill until the motor has completely stopped.
- Do not apply excess pressure to the bur as it may break or bend or become difficult to remove.

 **CAUTION**

- After the bur/drill is locked in place, pull and push the bur/drill to make sure the bur/drill is locked.
- Depress the push button until it is almost parallel to the head cap.
- If removing the bur/drill is difficult, grip the bur/drill with fine pliers while depressing the Push-button, and gently remove the bur/drill. This situation can be avoided by frequently mounting and removing the bur/drill.

6. Spray Nozzle

Three irrigation methods are available for X-SG20L / SG20; External, Internal or both simultaneously, depending on the instrument and operation procedure.

6-1 External Spray Nozzle

Connect the Irrigation Tube to the External Spray Nozzle firmly. (Fig. 5).

6-2 Internal Spray Nozzle

Only for a drill with internal irrigation system.

1) Attach the Nozzle Holder. (Fig.6)


Fig. 5


Fig. 6

2) Connect the Irrigation Tube to the Internal Spray Nozzle firmly. (Fig.7)


Fig. 7

3) Insert the end of the Internal Spray Nozzle into the hole at the top of the handpiece head, then attach the Internal Spray Nozzle into the slot on the Nozzle Holder. (Fig.8)


Fig. 8

6-3 Simultaneous External and Internal Irrigation

- 1) Connect one end of the Y-Connector to the External Spray Nozzle and the other end into the Internal Spray Nozzle. (Refer to "6-1, 6-2")
- 2) Connect the Irrigation Tube to the Y-Connector firmly.(Fig.9)


Fig. 9

7. Checking the Handpiece Before Each Use

Follow the check below before use. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer.

- 1) Check the Head Cap is firmly tightened.
- 2) Check the coolant water is flowing properly.
- 3) Insert the bur/drill into the handpiece. Then rotate the handpiece for about one minute with coolant water at the Max. rotation speed of the attached motor. During rotation, check for abnormalities such as abnormal rotation, vibration, noise.
- 4) After the handpiece rotation has completely stopped, touch the handpiece head to confirm the head is NOT heating abnormally (Fig. 10).


Fig. 10

CAUTION

- To avoid injury, keep your hands away from the bur during rotation.

8. Post-use Maintenance

After each patient maintain the product as follows.

WARNING

- Immediately after a surgery, the handpiece should be cleaned, lubricated and sterilized. Failure to properly maintain the handpiece can lead to contaminants remaining inside the handpiece. Unremoved blood may result in product failure due to the internal blood coagulation. Improper maintenance may also lead to overheating, causing burn injuries, product failure, or difficulty in attaching and detaching the bur/drill.

CAUTION

- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.

8-1 Preparation

- 1) Wear a eye protection, a mask and gloves to avoid infection.
- 2) Remove the bur.
- 3) (Only for X-SG20L / SG20) Immerse about 20mm of the handpiece head in water, and rotate the handpiece for about 15 seconds while switching the Forward/Reverse rotation direction (Fig. 11)
- 4) Disconnect the handpiece from the motor.


Fig. 11

8-2 Cleaning

■ Manual Cleaning (External Surface)

- 1) Remove dirt and debris on the handpiece surface using a brush under running water (at 38°C or less). Do not use a wire brush.
- 2) Wipe off the moisture on the surface using absorbent cloth.

CAUTION

- If water comes inside the handpiece, blow air into the rear of the handpiece using an air syringe and remove the remaining water.

■ Automatic Cleaning (External and Internal Surfaces)

Use a thermo-disinfector in accordance with ISO 15883 (EN ISO 15883).

Refer to the thermo-disinfector manual.

CAUTION

- After washing with thermo-disinfector and prior to lubrication, dry the product until all internal moisture is thoroughly removed. Thermo-disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product.

8-3 Cleaning (Optic)

Wipe clean the Glass Rod tip with an alcohol-immersed cotton swab. Remove all debris and oil.


Fig. 12

CAUTION

- Do not use a sharp tool to clean the Glass Rod. It could damage the glass and reduce the light transmission.

8-4 Cleaning (Spray Port) (X-SG20L / X-SG25L / SG20)

- 1) Internally clean the Internal Spray Nozzle and/or the External Spray Nozzle by using Cleaning Wire. (Fig. 13,14)
- 2) Connect the Irrigation Tube to each Spray Nozzle, and wash the inside with clean water. (Refer to "6. Spray Nozzle")


Fig. 13


Fig. 14

 **CAUTION**

- Do not forcibly insert the cleaning wire into the Spray Nozzle. Resultant port damage could cause the spray to be directed away from the bur/drill, and cause reduction of cooling efficiency.

8-5 Disinfection

■ Manual Disinfection (External Surface)

Wipe the surface with cloth soaked in disinfectant solution such as alcohol.

■ Automatic Disinfection (External and Internal Surfaces)

Use a thermo-disinfector in accordance with ISO 15883 (EN ISO 15883).

Refer to the thermo-disinfector manual.

8-6 Lubrication

■ NSK PANA SPRAY Plus / PANA SPRAY

Apply NSK PANA SPRAY Plus / PANA SPRAY every time after each use and/or before autoclaving.

- 1) Attach the Spray Nozzle into the Spray Port nozzle on the can.
- 2) Insert the Spray Nozzle into the rear of the handpiece. Hold the handpiece and spray 2-3 times for approximately 2-3 seconds each time. Apply lubricant until it expels from the handpiece head for at least 2 seconds (Fig. 15).


Fig. 15

⚠ CAUTION

- Hold the spray can upright.
- Be sure to hold the handpiece firmly to prevent slipping out of hand by the spray pressure when lubricating.
- Supply lubricant until it comes out of the handpiece head.
- If you wish to expel excessive oil from inside of the handpiece, rotate the handpiece for approx. 15 seconds without a bur. During the rotation, do not depress the push button inadvertently, especially when using absorbent cloth to prevent oil from scattering. Depressing the push button during rotation may cause a failure to attach a bur/drill.


Fig. 16


Fig. 17

NOTICE

- NSK recommends the use of "Spray Mist Absorber" (Y900084) to prevent oil mist expelling out of the handpiece head.

■ Chuck cleaning (X-SG93L / X-SG93)

Clean the Push Button chuck once a week.

- 1) Mount the Tip Nozzle into the spray can port.
- 2) Lubricate the chuck directly through the bur insertion hole (Fig. 18).
- 3) Lubricate handpiece by using NSK PANA SPRAY Plus / PANA SPRAY (Fig. 15) or NSK automatic handpiece cleaning and lubrication system.


Fig. 18

⚠ CAUTION

- If the chuck is not regularly cleaned the chuck grip may be weakened and the bur may be accidentally released during use.

8-7 Sterilization

Sterilize the product by autoclave sterilization. After each patient sterilize as noted below.

The following items can be autoclaved: Handpiece, Internal Spray Nozzle, Nozzle Holder

1) Insert into an autoclave pouch. Seal the pouch.

2) Autoclavable under the conditions below.

Exposure at 132°C for more than 15min.. Drying for more than 30 min.

3) The handpiece should remain in the autoclave pouch until required for use.

CAUTION

- Be sure to clean and lubricate the handpiece prior to autoclaving. If blood remains on the internal surface it can become clotted and cause product failure.
- Do not autoclave the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- To avoid product failure, do not use a sterilizer whose temperature exceeds 138°C during use including the dry cycle. The chamber temperature of some type of sterilizers may exceed 138°C. Contact the sterilizer manufacturer for detailed information.
- Keep the product in suitable atmospheric pressure, temperature, humidity, ventilation, and sunlight. The air should be free from dust, salt and sulphur.
- Do not touch the product immediately after autoclaving as it will be very hot and must remain in a sterile condition.
- Autoclave sterilization is recommended for the product. The validity of other sterilization methods is not confirmed.
- Only Handpiece, Internal Spray Nozzle and Nozzle Holder can be autoclaved. Do not autoclave other components (Y-Connector, Irrigation Tube, etc.).

NOTICE

- NSK recommends Class B sterilizers as stated in EN 13060.

9. Periodical Maintenance Checks

Perform periodical maintenance checks every three months, referring to the check sheet below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Head cap is loose	Check that the Head Cap is firmly tightened.
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.

10. Specifications

Model	X-SG20L	X-SG25L	X-SG93L	X-SG93	SG20
Max. Rotation Speed (Motor)	40,000min ⁻¹				
Max. Rotation Speed (Handpiece)	2,000min ⁻¹	40,000min ⁻¹	120,000min ⁻¹		2,000min ⁻¹
Gear Ratio	20:1 Reduction	1:1 Direct Drive	1:3 Increasing		20:1 Reduction
Bur/Drill Type	ISO 1797-1 (EN ISO 1797-1) Type1 Ø2.35mm Surgical Bur/Drill	ISO 1797-1 (EN ISO 1797-1) Type1 Ø2.35mm CA Bur	ISO 1797-1 (EN ISO 1797-1) Type3 Ø1.60mm FG Bur		ISO 1797-1 (EN ISO 1797-1) Type1 Ø2.35mm Surgical Bur/Drill
Chuck Length	11.6mm				
Max. Bur/Drill Length	36mm	22.5mm	25mm		36mm
Max. Working Part Diameter	Ø4.7mm	Ø4.0mm	Ø2.0mm		Ø4.7mm
Optic	Glass Rod	-	Glass Rod	-	
Water Spray Type	External, Internal*	External			External, Internal*
Use Environment	Temperature: 10 - 40°C, Humidity: 30 - 75% (No Condensation)				
Transportation and Store Environment	Temperature -10 - 50°C, Humidity: 10 - 85%, Atmospheric Pressure: 500 - 1,060hPa				

*Only for a drill with internal irrigation system.

11. Symbol

This product is Autoclavable up to Max.135°C.

This product can be washed via Thermo Disinfectant.

Conforms to CE European Directive of "Medical device directive 93/42/EEC."

Manufacturer.

Authorized representative in the European community.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

12. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

13. Spare Parts List

Model	Order Code	Compatible Product
Internal Spray Nozzle	10000324	X-SG20L / SG20
Nozzle Holder	20000396	X-SG20L
Nozzle Holder	20000357	SG20
Y-Connector	C823752	X-SG20L / SG20
Cleaning Wire	20000512	X-SG20L / X-SG25L / SG20
E Type Spray Nozzle	Z019090	-

14. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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2015.03.20 00 ©
NSK America Corp www.nskdental.com