SGL70M/SG70M

Please read this Operation Manual carefully before use, and file for future reference.

MADE IN JAPAN **ENGLISH** OM-E0820E 000

◆ Indications for use

The Surgic Pro⁺ / Surgic Pro is intended for use in dental oral surgery and dental implant. The main unit is designed to be used with a specific dental micromotor that drives dental handpieces fitted with appropriate tools to cut hard tissues in the mouth.

⚠ Cautions for handling and operation

- Read these safety cautions thoroughly before use and operate the product properly.
- These indicators are to allow you to use the product safely, prevent danger and harm to you and others. These are classified by degree of danger, damage and seriousness. All indicators concern safety, be sure to follow them.

◆ Classification of equipment

- Method of sterilization or disinfection recommended by the manufacturer :
-Refer to the operation manual for Surgic Pro series.

Classification	Degree of Danger or Damage and Seriousness
⚠ CAUTION	Explains an instruction where minor to medium injury or physical damage may occur.

1. Safety precautions prior to use

⚠ CAUTION

- This Micromotor is used only for the NSK Surgic Pro series. Do not connect this Micromotor to other Control Units.
- Do not attempt to disassemble the product nor temper with the mechanism.
- Do not disconnect the Motor Cord from the motor.
- Do not drop, hit, or subject to excessive shock.
- Do not wipe or clean with high acid water. Do not immerse in any solutions.
- This product is delivered in a non-sterile condition and must be autoclaved prior to use.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

2. Package Contents

No.	Description	Quantity
1	Optic Motor / Non-Optic Motor (With Motor Cord)*	1
2	Protection Plug	1
3	O-ring (Protection Plug)	2
4	O-ring (Motor)	2

*Either Motor is included.

3. Installation

Align the (●) mark on the Motor Cord with (●) mark on the Control Unit.

Push the Base of Motor Cord Plug into the socket until the locking joint 'clicks' into place.

To disconnect the plug; pull back the Lock Joint, pull further to disconnect the cord.


Fig. 1

4. Maintenance (Cleaning) and Sterilization

Refer to the operation manual for Surgic Pro series.

5. Specification

Model	SGL70M	SG70M
Type	E1023	E1025
Rotation Speed Range	200 – 40,000 min ⁻¹	
Input Voltage	DC24V	
Dimensions	Ø23.5xL103.3mm (without cord)	
Cord Length	2m	
Optic	White LED	-
Weight	Approx.280g	


	Temperature	Humidity	Atmospheric pressure
Use Environment	Between 0-40°C (32-104°F)	Between 10-85% RH	Between 700-1060 hPa
Transportation and Store Environment	Between -10-50°C (14-122°F)	Between 10-85% RH	Between 500-1060 hPa

No moisture condensation in the motor. Use at outside of these limits may cause malfunction.

6. Disposing Product

Consult with your dealer from whom you purchased the device regarding waste disposal.

7. Symbols

-  The EU directive 93/42/EEC was applied in the design and production of this medical device.
-  Caution, See Operation Manual
-  Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician
-  Manufacturer
-  This product can be sterilized in a steam sterilizer up to Max. 135°C.
-  The surface of the product can be used with a Thermo-Disinfectant prior to Sterilization.
-  GS1 DataMatrix for Unique Device Identifier.

8. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

Manufacturer

NAKANISHI INC.
www.nsk-inc.com

700 Shimohinata, Kanuma
Tochigi 322-8666, Japan

NSK America Corp
www.nskdental.com

1800 Global Parkway, Hoffman Estates
IL 60192, USA

Specifications are subject to change without notice.

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