The NSK logo consists of the letters "NSK" in a bold, italicized, blue sans-serif font. The letters are slightly slanted to the right.A rectangular icon with rounded corners containing the text "Rx Only" in a blue, sans-serif font. The icon has a light blue background.

Please read this Operation Manual carefully
before use, and file for future reference.

Surgical Handpiece

Ti-Max **X** X-SG20L
S-Max SG20

MADE IN JAPAN OM-C0560E 001

 CAUTION for X-SG20L / SG20

- This handpiece is intended for use only with surgery system with torque calibration function such as NSK Surgic Pro. Using this product with a different manufacturer's system may result in injuries or damage the handpiece.

1. User and Indications for Use

User: Qualified Professionals

Indications for Use: This medical device is for oral surgery and dental implant operation. This device is driven by an electronic micromotor for oral surgery and dental implant. The device is intended to transmit the rotation of the power source at different gear ratios, thereby moving instruments such as surgical burs or surgical drills to cut the maxilla/mandible during oral surgery and dental implant.

2. Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
 WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
 CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

WARNING

- Clean, disinfect, and lubricate the handpiece immediately (within 1 hour) after each treatment to remove residue. Failure to properly maintain the handpiece may cause infection, product failure, overheating leading to burn injuries. (Refer to "8. Post-use Maintenance")
- If blood infiltrates inside a handpiece, an automatic handpiece cleaning and lubrication system may not totally clean the internal handpiece components, which may lead to internal blood coagulation. Coagulated blood inside the handpiece may cause handpiece failure and overheating, causing burn injuries. To avoid the risk, NSK recommends to immediately lubricate those handpieces with PANA SPRAY Plus.
- Depressing the Push Button while the handpiece is in rotation may lead to overheating, causing burn injuries or product failure. Avoid the push button to contact with any oral tissue.
- Supply coolant water through the spray nozzle while using the handpiece. Other water supply methods than recommended may lead to overheating causing burn injuries, or product failure.
- Do not use the handpiece for polishing purposes. Polishing paste could enter the handpiece, causing the failure of the Push Button or burn injuries due to the overheating of the handpiece head.

 **WARNING**

- Keep any debris or other foreign materials away from inside the gear or the handpiece. Foreign materials remaining inside may lead to overheating, causing burn injuries or other accidents.
- Ball bearings are wear items and may require replacement. Be sure to inspect the handpiece before use (Refer to "7. Checking the Handpiece Before Each Use"). If abnormal vibration or noise are found, the bearings or other internal parts may be worn out or damaged. Damaged parts may lead to handpiece overheating causing burn injuries. If overheating is suspected, stop using the handpiece immediately and contact your Authorized NSK Dealer.

CAUTION

- When operating the product always consider the safety of the patient.
- Users are responsible for the operational control, maintenance and continual inspection of this product.
- Prior to clinical use, inspect the handpiece. Check for vibration, noise and overheating. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer. (Refer to "7. Checking the Handpiece Before Each Use")
- Do not disassemble or alter the handpiece except as recommended by NSK in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product. Deformation may cause the handpiece to fail during use.
- Operators and all others in the area must wear eye protection and a mask when operating this handpiece.
- Do not use bur/drills with problems listed below as the bur/drill may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur/drill.
 - Bur/drill which is cracked on the edge or axis.
 - Non-ISO(EN ISO) standard, or tampered bur/drill.
- Do not use bur/drills other than specified in "12. Specifications" as such burs may be accidentally released or break during rotation.

CAUTION

- Insert a sterilized bur/drill and always keep the bur/drill shank clean. When changing burs/drills, make sure that there is no dirt or debris in and around the bur/drill insertion hole of the handpiece. Debris in the chuck could cause infection, heating, poor bur concentricity or low chuck retention force.
- Always insert the all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur may occur.
- Do not exceed Max. Bur/Drill Length shown in "12. Specifications".
- Always follow the instructions provided by the bur/drill manufacturer.
- Do not exceed the bur/drill speed recommended by the bur manufacturer.
- Should the handpiece function abnormally during use, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- Do not use the following fluids to wipe, immerse or clean the product; strong/super acid water, strong acid/alkaline chemicals, chlorinecontaining solutions, solvents such as benzine or thinner.
- The handpiece is delivered in a non-sterile condition and perform steam sterilization prior to use.
- Perform regular function and maintenance checks. (Refer to "11. Periodical Maintenance Checks")

 **CAUTION**

- If the handpiece has not been used for a long period, rotate the handpiece and check for noise, vibration and overheating before use.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.
- The operation of the handpiece is permitted only on dental units which correspond to the standards IEC 60601-1 (EN 60601-1) and IEC 60601-1-2 (EN 60601-1-2).
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

NOTICE

- Repairs of this product are only to be performed by authorized service technicians according to NSK requirements. Contact your Authorized NSK Dealer if repairs are necessary.

3. Accessory List

No.	Part Name	Quantity
1	Internal Spray Nozzle	1
2	Nozzle Holder	1
3	Y-Connector	1
4	Cleaning Wire	1
5	E-Type Spray Nozzle	1

4. Connecting and Disconnecting the Handpiece from the Motor

4-1 Connecting

- 1) Insert the handpiece direct to the motor.(Fig. 1) When connecting an optic handpiece, screw the handpiece until it clicks and locks into position.
- 2) Confirm that the handpiece is firmly connected to the motor.


Fig. 1

4-2 Disconnecting

Hold the motor and the handpiece, then pull apart.

CAUTION

- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Connect ONLY to E type motors (ISO 3964 (EN ISO 3964)).
- Do not exceed the Max. Rotation Speed (Motor) shown in "12. Specifications".

CAUTION

X-SG20L / SG20

- This handpiece is intended for use only with surgery system with Torque Calibration System such as NSK Surgic Pro. Using this product with a different manufacture's may result in injuries or damage the handpiece.

SG20

- The handpiece is designed for connection to a NSK Non-optic motor or to a motor with a stopper pin. Connection to any other type of motor cannot lock the handpiece.

5. Mounting and Removing the Bur/Drill

5-1 Mounting the Bur/Drill

- 1) Insert the bur/drill until it is correctly seated in place.
- 2) Depress the Push Button and insert the bur into the chuck until the bur "notch" mechanism engages. Release the button.
- 3) Confirm the bur/drill is secure by gently pulling and pushing the bur/drill without depressing the Push Button.

5-2 Removing the Bur

Depress the Push Button to open the chuck (❶) and remove the bur/drill (❷).


Fig. 2


Fig. 3

 **WARNING**

- Always use a clean, sterile bur/drill. Be sure to wear sterile gloves before mounting the bur/drill.
- When changing burs/drills during treatment, remove the bur/drill and wipe around the bur/drill insertion hole and gloves with a dry lint free cloth until the insertion hole is visually clean, then mount a clean, sterile bur/drill. Residual blood or debris may enter the handpiece while changing burs, which could lead to higher risk of cross-contamination.

 **CAUTION**

- Do not use bur/drills with problems listed below as the bur/drill may break, seize up or disengage from the chuck.
 - Bent, deformed, worn, rusted, broken, deficient bur/drill.
 - Bur/drill which is cracked on the edge or axis.
 - Non-ISO(EN ISO) standard, or tampered bur/drill.
- Do not use bur/drills other than specified in "12. Specifications" as such burs may be accidentally released or break during rotation.

CAUTION

- Insert a sterilized bur/drill and always keep the bur/drill shank clean. When changing burs/drills, make sure that there is no dirt or debris in and around the bur/drill insertion hole of the handpiece. Debris in the chuck could cause infection, heating, poor bur concentricity or low chuck retention force.
- Always insert the bur/drill all the way into the chuck. If insertion is insufficient, premature failure of the bearings or accidental release of the bur may occur.
- Do not exceed Max. Bur/Drill Length shown in "12. Specifications".
- Always follow the instructions provided by the bur/drill manufacturer.
- Do not exceed the bur speed recommended by the bur/drill manufacturer.
- Do not mount or remove the bur/drill until the motor has completely stopped.
- Do not apply excess pressure to the bur as it may break or bend or become difficult to remove.
- After the bur/drill is locked in place, pull and push the bur/drill to make sure the bur/drill is locked.
- Depress the push button until it is almost parallel to the head cap.
- If removing the bur/drill is difficult, grip the bur/drill with fine pliers while depressing the Push-button, and gently remove the bur/drill. This situation can be avoided by frequently mounting and removing the bur/drill.

6. Spray Nozzle

Three irrigation methods are available; External, Internal or both simultaneously, depending on the instrument and operation procedure.

6-1 External Spray Nozzle

Connect the Irrigation Tube to the External Spray Nozzle firmly. (Fig. 4)

6-2 Internal Spray Nozzle

Only for a drill with internal irrigation system.

1) Attach the Nozzle Holder. (Fig. 5)


Fig. 4


Fig. 5

2) Connect the Irrigation Tube to the Internal Spray Nozzle firmly. (Fig. 6)


Fig. 6

3) Insert the end of the Internal Spray Nozzle into the hole at the top of the handpiece head, then attach the Internal Spray Nozzle into the slot on the Nozzle Holder. (Fig. 7)


Fig. 7

6-3 Simultaneous External and Internal Irrigation

- 1) Connect one end of the Y-Connector to the External Spray Nozzle and the other end into the Internal Spray Nozzle. (Refer to "6-1, 6-2")
- 2) Connect the Irrigation Tube to the Y-Connector firmly. (Fig. 8)


Fig. 8

7. Checking the Handpiece Before Each Use

Follow the steps below before use. If any abnormalities are found, stop using the handpiece immediately and contact your Authorized NSK Dealer.

- 1) Check the Head Cap is firmly tightened.
- 2) Check the coolant water is flowing properly.
- 3) Mount the bur/drill. (Refer to "5. Mounting and Removing the Bur/Drill")
- 4) Rotate the handpiece for about one minute with coolant water at the Max. rotation speed of the attached motor.
During rotation, check for abnormalities such as abnormal rotation, vibration, noise, water leakage.
- 5) After the handpiece rotation has completely stopped, touch the handpiece head to confirm the head is NOT heating abnormally. (Fig. 9)


Fig. 9

CAUTION

- To avoid injury, keep your hands away from the bur during rotation.

8. Post-use Maintenance

After each patient maintain the product as follows. Post-use maintenance includes cleaning, lubricating, and sterilization after each patient use.

WARNING

- To ensure cleaning and sterilization efficacy, use only the following procedures for reprocessing.
- Follow any additional local directives, standards, and guidelines for cleaning, disinfection, and sterilization.
- Clean, disinfect, and lubricate the handpiece immediately (within 1 hour) after each treatment to remove residue. Failure to properly maintain the handpiece may cause infection, product failure, overheating leading to burn injuries.

8-1 Cleaning at point-of use

CAUTION

- Do not use the following liquids to wipe, immerse, or clean the product: strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.
- Do not immerse NSK instruments in disinfectant solution or clean in ultrasonic devices.

1) Always wear protective gloves, a mask, and protective goggles for safety purposes and to minimize the risk of infection. (Fig. 10)


Fig. 10

2) Remove the bur/drill. (Fig. 11)


Fig. 11

3) Operate the handpiece while running clean water at maximum flow rate into the irrigation tube with the irrigation nozzle, at the chair-side for at least 20 seconds to purge fluids in the handpiece (Fig. 12). After the handpiece rotation has completely stopped, touch the handpiece head to confirm the head is NO heating abnormally. (Water temperature 38°C or less, drinking water level of water quality recommended)


Fig. 12

⚠ CAUTION

- If the head is heating abnormally, stop using the handpiece immediately and contact your Authorized NSK Dealer.

- 4) Remove the irrigation tube. Immerse about 20mm of the handpiece head in water then rotate the handpiece in forward and reverse rotation at maximum rotation speed for about 15 seconds. (Fig. 13)


Fig. 13

- 5) Wipe the exterior of the handpiece clean using a cloth moistened with disinfectant or ethanol with a concentration of 60-90%. (Fig. 14)

When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties.

The following disinfectants can be used in the United States and Canada.

CaviCide[®], CaviWipes[®] (manufactured by Metrex)


Fig. 14

- 6) Remove the handpiece from the motor. (Fig. 15)

Then, remove the internal spray nozzle and nozzle holder from the handpiece.

- 7) Carry the handpiece and the internal spray nozzle and nozzle holder to the decontamination area.


Fig. 15

8-2 Cleaning, Disinfecting, and Drying

Perform cleaning, disinfecting, and drying using either the manual or automatic method.

■ Manual Method (Cleaning, Disinfecting, and Drying the Exterior)

Wash the handpiece and irrigation nozzle under running water following the procedures below.

(Water condition: $\leq 38^{\circ}\text{C}$, $\geq 3.5\text{L}/\text{min}$; water should be of the same quality as drinking water)

- 1) Clean the external surfaces of the handpiece with a soft bristled toothbrush for 15 seconds or more. (Fig. 16)


Fig. 16

- 2) Clean inside the bur/drill insertion hole for 15 seconds or more, using an interdental brush (wire thickness 0.7mm) that conforms to size 4 in ISO 16409. (Fig. 17)


Fig. 17

- 3) Clean the external surfaces of the internal spray nozzle and nozzle holder with a soft bristled toothbrush for each 15 seconds or more. (Fig. 18)


Fig. 18

4) After using the cleaning wire, clean the wire by pulling both ends of it with fingers while running water all around the wire for 15 seconds or more (Fig. 19).

Under appropriate lighting (natural office lighting, 500 lx or higher), inspect the handpiece, internal spray nozzle, nozzle holder and cleaning wire for blood or other debris. If any visible debris remains, repeat the process until they are visually clean.


Fig. 19

5) Wipe the exterior of the handpiece, internal spray nozzle, nozzle holder and cleaning wire clean using a cloth moistened with disinfectant or ethanol with a concentration of 60-90%. (Fig. 20)


Fig. 20

When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties.

The following disinfectants can be used in the United States and Canada. CaviCide[®], CaviWipes[®] (manufactured by Metrex)

After cleaning and disinfection, if the residual moisture is present, wipe off with dry cloth or blow it off with compressed air ($\leq 0.35\text{MPa}$) until there is no moisture in the interior and exterior. (Only for the cleaning wire, wipe off with dry cloth.)

When blowing off with compressed air, cover the handpiece with cloth to prevent scattering of water. Proceed to "8-3 Cleaning and Lubrication".

Automatic Method (Cleaning, Disinfecting, and Drying the Exterior and Interior)

 NSK handpieces with this symbol are compatible with medical thermo-disinfectors.

Due to the variation in cleaning/disinfection cycles and cleaning agents available from the various manufacturers, NSK has validated the Getinge 46 Series Washer Disinfector Cycle P7-Anesthesia for automatic processing/reprocessing of the product (handpiece, internal spray nozzle, nozzle holder and cleaning wire).

Only validated procedures must be used for cleaning, disinfection and sterilization. When using reprocessing procedures that are different from those described in this manual, those procedures must be validated by the respective practice or hospital using a thermo-disinfector that conforms to ISO 15883-1.

NOTICE

- All water is "building supply" (at least drinking water quality) unless otherwise indicated.

- 1) Prepare and pre-disinfect the device as described above in step "8-1 Cleaning at point-of use".
- 2) Clean inside the bur/drill insertion hole for 15 seconds or more under running water, using an interdental brush (wire thickness 0.7mm) that conforms to size 4 in ISO 16409. (Fig. 21) (Water condition: $\leq 38^{\circ}\text{C}$, $\geq 3.5\text{L/min}$; water should be of the same quality as drinking water)
- 3) Wipe off a dirt on a part held by a handpiece holder using the cloth moistened with 60 - 90% ethanol or disinfectant.


Fig. 21

When using disinfectant, follow the instructions given by the manufacturer of the disinfectant.

Use a state-sanctioned disinfectant with proven bactericidal, fungicidal, and virucidal properties.

The following disinfectants can be used in the United States and Canada.

CaviCide[®], CaviWipes[®] (manufactured by Metrex)

- 4) Perform main wash under the following conditions. Use a Getinge 46 Series thermo-disinfector (or equivalent device):
 - Pre-Wash 1: Cold water (<110°F (43°C)) for 2 min
 - Pre-Wash 2: None - Not Applicable
 - Enzyme Wash: Hot tap water and 60ml of heated (95°F (35°C)) detergent (Getinge Clean Renuzyme) for 3 min at 120°F (49°C)
 - Detergent Wash: None - Not Applicable
 - Rinse 1: Hot tap water for 2 min
 - Rinse 2: Hot tap water for 2min
 - RO Final Rinse: Hot tap water for 32 min at 167°F (75°C)
 - Drying: 20 minUpon removal from the thermo-disinfector, if the residual moisture is present, wipe off with dry cloth or blow it off with compressed air ($\leq 0.35\text{MPa}$) until there is no moisture in the interior and exterior. (Only for the cleaning wire, wipe off with dry cloth.)
- 5) After cleaning, disinfecting, and drying the product, confirm that it is clean under appropriate lighting (500 lx or higher). If dirt still remains, repeat the procedures.
- 6) Proceed to "8-3 Cleaning and Lubrication".

CAUTION

- Handpieces must be removed from the thermo-disinfector immediately (within 1 hour) after the cleaning, disinfecting and drying cycle is complete to prevent corrosion.
- After washing with thermo-disinfector, prior to lubrication, dry the product until all internal moisture is completely removed. Thermo-disinfector moisture remaining inside the product could reduce the effect of lubrication and could cause corrosion inside of the product.

8-3 Cleaning and Lubrication

Perform cleaning and lubrication using the manual method.

Manual Method (Cleaning and Lubricating the Interior)

- 1) Attach an E-Type Spray Nozzle to the nozzle of the PANA SPRAY Plus.
- 2) Shake the can 3 or 4 times, and insert the E-Type spray nozzle into the rear of the handpiece while holding the head of the handpiece with a cloth etc..
- 3) Hold the handpiece, and spray for 2 or 3 seconds until oil comes out of the handpiece head. (Fig. 22)


Fig. 22

Repeat lubrication until foreign material stops coming out of the tip.

⚠ CAUTION

- Do not use sprays other than NSK PANA SPRAY Plus. Using sprays other than PANA SPRAY Plus may cause the handpiece to overheat.
- Hold the spray can upright.
- Firmly hold the handpiece to prevent it slipping when spray pressure is applied.
- Spray lubricant until it expels from the handpiece head.
- If you wish to purge excessive oil from inside the handpiece, rotate the handpiece for approx 15 seconds without a bur/drill. During the rotation, do not depress the push button, especially when using absorbent cloth to prevent oil from scattering. Depressing the push button during rotation may damage the chuck mechanism.


Fig. 23


Fig. 24

NOTICE

- NSK recommends the use of "Spray Mist Absorber" (Y900084) to prevent oil mist expelling out of the handpiece head.

8-4 Packaging, Sterilizing, and Drying

- 1) Insert the handpiece, internal spray nozzle, nozzle holder and cleaning wire into an FDA-approved sterilization pouch that conforms to ISO 11607-1, and seal the pouch.
- 2) After packaging in a sterile pouch, sterilize using the following validated cycles. Use an FDA-cleared steam sterilizer to perform sterilizations.

Type	Gravity Displacement	Pre-Vacuum (Dynamic Air Removal)
Temperature	132°C	132°C
Full Cycle Time	15 min	4 min
Drying Time	30 min	30 min

 **CAUTION**

- Use an FDA-cleared steam sterilizer to perform sterilizations.
- Follow local rules, regulations, and guidelines regarding the reprocessing of devices.
- Do not touch the product immediately after steam sterilization as it will be very hot and must remain in a sterile condition.
- Do not perform steam sterilization on the product with other instruments even when it is in a pouch. This is to prevent possible discoloration and damage to the product from chemical residue on other instruments.
- Clean and lubricate the handpiece prior to sterilization. If blood remains on the internal surface it can become clotted and cause product failure.
- Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product.
- Be sure to use sterilizers that can perform sterilization up to 135°C. In some sterilizers, the chamber temperature may exceed 135°C. Do not use these sterilizers as failure of the handpiece could occur. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Steam sterilization is recommended for the product. The validity of other sterilization methods (such as plasma sterilization or EOG sterilization) is not confirmed.

8-5 Storage

Store the product in a dry, clean location.

CAUTION

- After the sterilization and drying cycles are complete, remove the handpiece immediately from the sterilizer to store it.
- Store the product in a well ventilated place out of direct sunlight and within the range of temperature, humidity and pressure specified in "12. Specifications".
- Sterilization is not guaranteed after the sterilization retention period specified by the manufacturer and seller of the sterilization pouch has elapsed. If the sterilization retention period has elapsed, perform sterilization again with a new sterilization pouch.

9. Cleaning of Optic Illumination Points

When debris is attached to optic illumination points (Glass Rod), wipe clean all the optic illumination points using an alcohol-immersed cotton swab. Remove all debris.

CAUTION

- Do not use a sharp tool to clean the Glass Rod. It could damage the glass and reduce the light transmission.

10. Cleaning of Spray Nozzle

When spray nozzle are clogged, or spray does not exit evenly from the nozzle, internally clean the internal spray nozzle and/or the external spray nozzle by using cleaning wire. Connect the irrigation tube to each spray nozzle, and wash the inside with clean water.


Fig. 26


Fig. 25


Fig. 27

After using the cleaning wire, reprocess the cleaning wire following "8. Post-use Maintenance".

 **CAUTION**

- Do not forcibly insert the cleaning wire into the Spray Nozzle. Resultant port damage could cause the spray to be directed away from the bur/drill, and cause reduction of cooling efficiency.

11. Periodical Maintenance Checks

Perform periodical maintenance checks every three months, referring to the check sheet below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Head cap is loose	Check that the Head Cap is firmly tightened.
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.

12. Specifications

Model	X-SG20L	SG20
Max. Rotation Speed (Motor)	40,000min ⁻¹	
Max. Rotation Speed (Handpiece)	2,000min ⁻¹	
Gear Ratio	20:1 Reduction	
Bur/Drill Type	ISO 1797-1 (EN ISO 1797-1) Type1 Ø2.35mm Surgical Bur/Drill	
Chuck Length	11.6mm	
Max. Bur/Drill Length	36mm	
Max. Working Part Diameter	Ø4.7mm	
Max. Torque	80N·cm	
Optic	Glass Rod	-
Water Spray Type	External, Internal*	
Use Environment	Temperature: 10 - 35°C, Humidity: 30 - 75% (No Condensation), Pressure: 700 - 1,060hPa	
Transportation and Store Environment	Temperature -10 - 50°C, Humidity: 10 - 85%, Pressure: 500 - 1,060hPa	

*Only for a drill with internal irrigation system.

13. Symbol

This product can be sterilized in a steam sterilizer up to Max.135°C.

This product can be washed via Thermo Disinfector.

Conforms to CE European Directive of "Medical device directive 93/42/EEC."

Manufacturer.

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

GS1 DataMatrix for Unique Device Identifier.

14. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

15. Spare Parts List

Model	Order Code	Compatible Product
Internal Spray Nozzle	10000324	X-SG20L / SG20
Nozzle Holder	20000396	X-SG20L
Nozzle Holder	20000357	SG20
Y-Connector	C823752	X-SG20L / SG20
Cleaning Wire	20000512	X-SG20L / SG20
E-Type Spray Nozzle	Z019090	-

16. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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NSK America Corp www.nskdental.com