OPERATION MANUAL

Rx Only

MADE IN JAPAN

OM-DE1011EN 001

Ultrasonic Handpiece

VA2-HP / VA2-LUX-HP / VA2-LUX-HP Ti

1 Preface

In this operation manual, VA2-LUX-HP Ti, VA2-LUX-HP, VA2-HP are called “ultrasonic handpiece”.

1-1 User and indications for use

-User

Dentists, oral surgeons, and hygienists

-Indications for use

Scaling, Perio, Implant Maintenance, Endodontic, Retrograde Endo, Restorative (for Minimal Intervention/Finishing/Trimming/Polishing/Caries of Dentin), Prosthetics (Condensation/Loosening/Plugging).

1-2 Related manuals

For proper use with the combined items, refer to the following operation manuals:

Operation manual Title	Operation manual No.
Varios 170 OPERATION MANUAL	OM-DE1010EN
Varios Tips OPERATION MANUAL	OM-DZ0919EN
CLA-10 OPERATION MANUAL	OM-DE0982MA
POWER GUIDE	OM-DE1247EN

1-3 Precautions for handling and operation

Read these cautions carefully and use only as intended or instructed.

Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device.

Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
⚠ CONTRA-INDICATIONS	Specific medical reason for not using the device : i.e. health condition of the patient.
⚠ WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
⚠ CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

⚠ CONTRA-INDICATIONS

- This handpiece is not intended for use on newborn (neonate) and infant (<2 years old) patient populations.

WARNING

- Patients (including infants) and users with external/implantable pacemakers and implantable cardioverter defibrillators. It may affect the operation of external/implantable pacemakers and implantable cardioverter defibrillators.
- Be sure to always carefully read the operation manual for the Ultrasonic System before use to fully understand the usage, performance, and functions of the instrument, and keep it for future reference.
- Do not use in rooms with a risk of explosion or near flammable substances. Also, do not use on or near patients under flammable anesthesia (nitrous oxide).
- Before unpacking this handpiece, make sure that the packaging has not been opened, tampered with, or damaged. If it has been opened or damaged, there may be a risk of malfunction or contamination. Stop use and contact an authorized NSK dealer.
- Always use irrigation in normal use (except for dry tips, which do not require coolant water). If there is insufficient irrigation, the product will overheat and may injure the tooth surface.
- Do not aim the light of the ultrasonic handpiece or handpiece cord at the eye of the patient, operator, or assistant. Doing so may damage eyesight.

⚠ WARNING

- Report all serious events, related to the product immediately to the manufacturer and the responsible national authority.
- Qualified professionals must use this device by fully complying with their respective country's regulations, accident prevention measures, and strictly follow these instructions for use. Users are responsible for operational control, maintenance, and inspection.
- When operating the product, always consider the safety of the patient.
- The user shall be responsible for any judgment that relates to the application of this handpiece to a patient.
- Should the product function abnormally, stop using the product immediately and contact an authorized NSK dealer.

CAUTION

- The handpiece surface may reach the following temperature when using a dry tip continuously for 23 minutes in a room temperature of 26 °C (78.8 °F). If the product becomes hot, stop use and leave it for approximately 10 minutes until it cools down.

Tip : 29 °C (84.2 °F)

Titanium handpiece (VA2-LUX-HP Ti) : 41.2 °C (106.2 °F)

Resin handpiece (VA2-LUX-HP, VA2-HP) : 50.2 °C (122.4 °F)

- Surface temperature of tip shall be more than 50 °C (122 °F) without using a coolant water. To avoid this event, be sure to use a coolant water.
- Use or store this handpiece in a location that meets the requirements for the use environment, storage and does not produce condensation. Otherwise, it may cause malfunction.

CAUTION

- The use of non-NSK tips may cause the following faults or unexpected accidents:
 - Inappropriate vibration caused by incompatible threads.
 - Accidental ingestion or aspiration of a broken tip by the patient.
 - Damage or premature wear of the ultrasonic handpiece threads.
- The glass rod of the ultrasonic handpiece can deteriorate as a result of a long term use over the lifetime of the device.
- Hold the connector part when pulling out the handpiece cord. Pulling the cord may cause a broken wire or damage.
- If water droplets, chemicals, solvents, or liquid disinfectant come into contact with this handpiece (including each part and each connection part), wipe them off immediately. Failure to do so may result in discoloration or deformation.
- Clean, disinfect, and sterilize it following the instructions in this manual before use.

CAUTION

- The effectiveness of sterilization methods other than autoclave sterilization is not verified.
- The ultrasonic handpiece must only be used by a qualified personnel in accordance with the intended purpose and the directions for use stated in the operation manual of the ultrasonic system of NSK and this operation manual. Do not use it for any other purpose.
- The ultrasonic handpiece is consumable. To ensure proper function and safe use of the medical device, perform periodical maintenance checks in accordance with the operation manual of the ultrasonic system.
- This handpiece is for indoor use only.
- Do not drop, hit, or subject to excessive shock. Glass rod and ultrasonic handpiece may be damaged.
- For information related to the EMC of this handpiece, refer to the operation manual of the ultrasonic system used.
- Push the handpiece securely and fully onto the handpiece cord when connecting.
- Before connecting the products, make sure that there is no dirt or foreign object adhered to each connection part. If there is any dirt or foreign object, clean it off.
- Read this operation manual before use to fully understand the product functions and file for future reference.

CAUTION

- Be sure to perform inspection by operating the product with each connecting part mounted before use.
- If any abnormal condition is detected, stop using the product and contact an authorized NSK dealer.
- If the product has not been used for a long period, check that it is functioning correctly and safely before using on a patient.
- Do not bend the handpiece cord at excessive angles. This may result in a decrease in the amount of irrigation.
- Immediately after using on a patient, clean, disinfect, and sterilize this handpiece with the correct method.
- Do not attempt to disassemble the product or tamper with the mechanism except as recommended by NSK in this manual.
- To avoid clinical downtime, it is recommended that a spare be kept on hand in case of a breakdown during operation.
- Be sure to conduct periodical maintenance check on the device and parts.
- This handpiece using for NSK ultrasonic system can be safely used in the treatment of children, pregnant or breastfeeding women, or other patient groups considered particularly vulnerable to such substances.
- Check that clean water is sprayed adequately before use.

CAUTION

- No special training is required for this device.
- U.S. Federal law restricts this device to sale by or on the order of a licensed physician.

1-4 Symbol

EN

FR

Use Environment

Move in the direction indicated.

Move fully to the stop, in the direction indicated.

Prohibited / Process NOT permitted / Wrong

Do the step / Process permitted / Correct

If...(Hypothetical situation)

Listen for the "click" sound.

Wear protective gloves

Wear a face mask

Wear eye protection

Temperature

The instrument can be used with a washer-disinfector for thermal disinfection.

This product can be sterilized in a steam sterilizer at 135 °C (275 °F).

Lamp / lighting / illumination

Drinking water

Visual inspection

Hand operated

Caution, hot surface

Consult operation instructions

Symbol for Caution.

Specifications

Manufacturer

Importer

Medical Device

Rx Only Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

Device Name for Health Canada

Catalog number (Order code)

Serial No.

Date of manufacture

Unique Device Identifier

GS1 DataMatrix for
Unique Device Identifier.

Contains hazardous
substances

Transportation and
Storage Environment

Temperature limitation

Humidity limitation

Atmospheric pressure
limitation

Keep away from
sunlight

Keep dry

fragile / handle with
care

This is the correct upright
position of the distribution
packages for transport
and/or storage

Auto Cleaning

NSK original surface
treatment

1–5 Disposing of products

When your medical device has reached the end of its service life, it must be disposed of via a specialist collection, removal, recycling or destruction channel. Prior to do so, in order to avoid the health risks to operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, you need to confirm the equipment is decontaminated.

Therefore contact an authorized NSK dealer.

1–6 Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model. Contact an authorized NSK dealer if repairs are necessary.

2 Product Description

2-1 Package contents

①

②

③

No.	Product Name	REF	Quantity
①	VA2-LUX-HP Ti	EA1203X050	1
②	VA2-LUX-HP	E351050	1
③	VA2-HP	E350050	1
-	Operation Manual	-	1

NOTICE

- One of items ① to ③ is included.

2-2 Applicable products

Product Name
Varios tip
Varios 170, Varios 970 (VA2-LUX-HP, VA2-HP, VA2-LUX-HP Ti)
Varios 370 (VA2-LUX-HP, VA2-HP)
Varios Combi Pro2 (VA2-LUX-HP Ti)

2-3 Part name

① VA2-LUX-HP Ti

② VA2-LUX-HP

③ VA2-HP

Handpiece

3 Preparation for Use

3-1 Mounting and Removing the handpiece

⚠ CAUTION

- Before connecting, make sure that dust or moisture does not adhere to the electrical contacts of the handpiece cord.
- Before connecting the handpiece to the handpiece cord, face the connector side of the handpiece towards a bright area and check the tip of the handpiece (for those with light). It is normal if the light-emitting part at the tip appears uniformly bright when illuminated. If part or all of it appears dark, the light-emitting part may be cracked or broken.
- Do not connect handpieces other than those listed in “2-2 Applicable products”.

Align the Dots on the handpiece and the handpiece cord. Push handpiece into connector.

To remove the handpiece, grip the handpiece and handpiece cord and pull to part handpiece and cord. (Figure 1)


Figure 1

⚠ WARNING

- To avoid an electric shock, do not touch the handpiece backend electrical contacts.

⚠ CAUTION

- Always confirm that the handpiece is correctly seated and locked into place.
- When removing the handpiece, always remove the tip first. There is a risk of injuring your hand with the tip.
- When attaching the handpiece to the handpiece cord, insert it firmly until it is fully seated. Attach and remove the handpiece and handpiece cord plug in parallel. Failure to do so may damage the inner surface of the handpiece cord plug and cause water leakage.
- When removing the handpiece from the handpiece cord, hold the handpiece cord plug to remove it. Pulling the cord may cause it to break. Do not twist the handpiece and cord to remove them. (Figure 2)


Figure 2

3-2 Using a barrier sleeve for handpiece cord

Use a disposable barrier sleeve that meets the U.S. and Canada regulations for the handpiece and handpiece cord, in order to prevent cross contamination between each patient. Install the disposable barrier sleeve to the handpiece cord (including the end of handpiece), then fix the barrier sleeve with a tape. (Figure 3)


Figure 3

⚠ CAUTION

- Do not wrap any tape or labels around the exterior of the handpiece other than tape for securing barrier film.
- During use or reprocessing, parts of the tape or labels may enter the handpiece, causing water blockages or vibration issues.
- If tape for securing barrier film is used, ensure it is completely removed after use.

3-3 Mounting and Removing the tip

- 1) Ensure that there is no debris on the threads of the handpiece and tip. If debris is present, remove it by blowing air or using another device. (Figure 4)


Figure 4

- 2) Turn the tip lightly by hand and install it. (Figure 5)


Figure 5

- 3) Insert the tip through the hole in the tip wrench.
(Figure 6)

Figure 6

- 4) Turn the tip wrench in the tightening direction indicated by the arrow until you hear a “click” sound, indicating it is securely tightened. (Figure 7)


Figure 7

- 5) To remove the tip, follow steps 3) to 4) above. Turn the tip wrench in the direction of the arrow to remove the tip. (Figure 8)


Figure 8

⚠ WARNING

- When replacing the tip, be careful if the tip protrudes from the tip wrench, because there is a risk of injury. Check the tip's operation manual for the appropriate tip wrench.
- Always use the NSK's tips with this product.
- Using tips from other manufacturers may cause the following issues or unexpected accidents:
 - Poor vibration due to forced installation of incompatible threads
 - Patient ingestion or aspiration due to tip breakage
 - Damage or premature wear of the handpiece threads
- Any issues or accidents caused by using tips from other manufacturers, including damage to this product, are not covered by the warranty.

CAUTION

- When tightening the tip, be careful not to twist the handpiece cord.
- Do not use tips that are damaged, deformed, or rusty, because they may break during use.
- If you turn the tip wrench at an angle, the threads may bend, the tip may not tighten properly, and vibration may be insufficient.
- If debris remains on the threads of the tip or handpiece during installation, it may cause abnormal noise or insufficient vibration. Clean the threads immediately if debris is present.
- Do not sharpen, bend, or change the angle of the tip, because it may break, produce abnormal noise, or result in insufficient vibration during use.
- Tighten the tip using the tip wrench or the E TIP WRENCH. If the tightening is weak, it may cause insufficient vibration or abnormal noise, or the tip may come off during operation, leading to accidents (refer to the tip's operation manual for the use of the E TIP WRENCH).
- When attaching the tip, wear clean gloves and use sterilized tips, handpieces, and tip wrenches to avoid injury.
- If the tip loosens, produces abnormal noise, or has weak vibration even after tightening with the tip wrench, the wrench may be worn out. The tip wrench is a consumable item and should be replaced approximately once a year.

CAUTION

- The G16 Tip can only be connected to the Varios2 handpiece (VA2-LUX-HP, VA2-HP, VA2-LUX-HP Ti). Using it with other handpieces may cause breakage.

4 Check before Operation

Before use, operate the device outside the patient's mouth to ensure there are no abnormalities such as vibration, noise, or heat generation, and check the following:

- The light turns on (for handpieces with lights).
- Water is properly sprayed from the tip.
- The tip vibrates properly.
- There are no looseness, abnormal vibrations, noises, or heat generation in the tip.

If you notice any abnormalities during inspection or use, stop using the device and contact an authorized NSK dealer.

5 Post-use Maintenance

Maintenance is performed after each patient treatment. The reprocessing instructions provided below are created based on the requirements of ISO 17664. The effectiveness of reprocessing is confirmed using the following procedures.

WARNING

- The product lifetime and the ability of the product to operate correctly are mainly determined by mechanical stresses during use and the chemical influences due to processing. If you detect any abnormality such as center heat generation, vibration or noise, or if the exterior markings (serial number or model number) have become unreadable before or during use, stop using the product immediately and ask an authorized NSK dealer for a periodic inspection.
- Disappearance of exterior markings is not covered under warranty for repair.
- When the product has been reprocessed 250 times or has been in use for one year, it is recommended that a periodic inspection be performed by an authorized NSK dealer.
- Failure to properly maintain the products may cause infection, product failure, or overheating and injuries. Be sure to clean and disinfect the products immediately (within 30 minutes) after use to remove residue.

WARNING

- Do not sterilize with ultraviolet rays. There is a risk of discoloration. The effectiveness of sterilization methods other than autoclave sterilization has not been confirmed.

CAUTION

- Wear eye protection, protective mask, and gloves for safety and infection prevention before starting the maintenance.
- Always obey the local directives, standards, guidelines, etc., as to the cleaning, disinfection, and sterilization of the product.
- Do not use solvents such as electrolyzed-oxidizing water (strong acidic water or super acidic water), strong acidic or alkaline chemical, chlorine-containing solvent, benzine and thinner when cleaning, immersing, or wiping the product.
- Ultrasonic handpiece is delivered in a non-sterile condition and requires reprocessing prior to first use and after each patient.

CAUTION

- Clean the ultrasonic handpiece prior to disinfection. If blood remains on the internal surface, it can become clotted and cause product failure.
- Please use disinfectant ethanol or disinfectant isopropanol (hereinafter referred to as disinfectant alcohol) that does not contain additives (such as quaternary ammonium salts) for cleaning and disinfection. Using other disinfectants may cause discoloration or cracking.
- Use appropriate disinfectant alcohol in accordance with the regulations and guidelines of the region where the product is used.
- Do not heat or cool the product too quickly. A rapid change in temperature could cause damage to the handpiece.
- Do not clean an ultrasonic handpiece with an ultrasonic cleaner.
- Use clean and dry compressed air for drying.
- This product requires reprocessing after each patient treatment.
- For reprocessing of the tips and tip wrench, please refer to the operation manual of the tips being used.

5-1 Preparation for maintenance

Items to prepare

- Eye protection
- Mask
- Gloves
- Cloth dampened with disinfectant alcohol
- FDA-approved sterilization pouch that conforms to ISO 11607-1

- 1 Wear eye protection, a mask, and gloves to prevent infection. (Figure 9)

Figure 9

- 2 Operate the tip at the Dental Delivery System-side for at least 20 seconds to purge fluids in the handpiece and tip.

- ③ Remove the tip.
To remove the tip, turn it counterclockwise with the tip wrench. (Figure 10)


Figure 10

For reprocessing of the tip, refer to the operation manual of the tip.

NOTICE

- If you are using Varios 970 or Varios Combi Pro2, perform auto-cleaning according to the Varios 970 or Varios Combi Pro2 operation manual. After auto-cleaning, proceed to the next step. If you are using other models, there is no auto-cleaning function. Proceed to the next step.

- ④ Remove the barrier sleeve and wipe the exterior and top of the handpiece with a damp, lint-free cloth for 10 seconds.
- ⑤ Remove the handpiece from the handpiece cord.

5-2 Pre-cleaning

- 1 Wash the surface and the top of the handpiece under running water using a brush (not metal) at a temperature below 38 °C (100.4 °F), with a flow rate of at least 3.5 L/min, for a total of 30 seconds.


Figure 11

- 2 Wash the gap between the handpiece's transducer threads and the top of the handpiece sheath under running water at a temperature below 38 °C (100.4 °F), with a flow rate of at least 3.5 L/min, for 40 seconds, using an interdental brush (wire thickness 0.7 mm) that conforms to size 4 in ISO 16409.


Figure 12

- ③ For VA2-LUX-HP / VA2-HP, rinse with running water for 10 seconds from the front side of the handpiece. Rinse with running water for 10 seconds from the connector side.

Then, hold the handpiece top downward and shake the handpiece twice to remove the moisture inside.

For VA2-LUX-HP Ti, rinse with running water for 10 seconds from the connector side.

Then, hold the handpiece top downward and shake the handpiece twice to remove the moisture inside.

In case of VA2-LUX-HP,
VA2-HP


Figure 13

In case of VA2-LUX-HP Ti


Figure 14

- ④ Wipe off any residual water on the exterior of the handpiece with a lint-free cloth.


Figure 15

- ⑤ Wipe the exterior of the handpiece with a lint-free cloth moistened with disinfectant alcohol.


Figure 16

- ⑥ Under appropriate lighting (500 lx or more), check that there is no dirt. If there is dirt, repeat from step 1 until the dirt is removed.

Proceed to sterilization. ⇒ 5-6 Packaging, Sterilizing, and Drying

Figure 17

5-3 Manual cleaning

- ① Prepare the cleaning solution in the cleaning tank according to the following conditions.
Prepare the cleaning solution according to the instructions of the detergent manufacturer.
Water quality: Distilled Water, Detergent: Medical enzyme-based neutral detergent (pH: 7.0 - 8.0)

- 2** Insert the handpiece into the cleaning solution. Then, angle the top of the handpiece upward and immerse it in the cleaning solution. Keep it submerged until all the air inside the handpiece is expelled.


Figure 18

- 3** Remove the handpiece from the cleaning solution and tilt the top downward to drain the cleaning solution.


Figure 19

- ④ Repeat step 2 and step 3 three times.


Figure 20

- ⑤ Perform step 2. Then, immerse the handpiece for 10 minutes.


Figure 21

- ⑥ Remove the handpiece from the cleaning solution and tilt the top downward to drain the cleaning solution. Repeat step 2 and step 3 three times.


Figure 22

- 7 For VA2-LUX-HP / VA2-HP, rinse with running water for 10 seconds from the front side of the handpiece. Rinse with running water for 10 seconds from the connector side.

Then, hold the handpiece top downward and shake the handpiece twice to remove the moisture inside.

For VA2-LUX-HP Ti, rinse with running water for 10 seconds from the connector side.

Then, hold the handpiece top downward and shake the handpiece twice to remove the moisture inside.


Figure 23

Figure 24

- 8 Wipe off any residual water on the exterior of the handpiece with a lint-free cloth.


Figure 25

- 9 Wipe the exterior of the handpiece with a lint-free cloth moistened with disinfectant alcohol.


Figure 26

- 10 Under appropriate lighting (500 lx or more), check that there is no dirt. If there is dirt, repeat from step 1 until the dirt is removed.

Proceed to sterilization. ⇒ 5-6 Packaging, Sterilizing, and Drying

Figure 27

5-4 Automatic processing/ reprocessing

- 1 Wipe the exterior of the handpiece connector with a lint-free cloth moistened with disinfectant alcohol.


Figure 28

- ② Attach the cleaning adapter (CLA-10: REF ZA11900102) to the washer-disinfector, and then attach the handpiece. The types of washer-disinfectors compatible with the CLA-10 and the gaskets for the washer-disinfectors are listed below.


Figure 29

Reference for washer-disinfector gasket *1

Manufacturer	REF	Connection Ø	Color
Miele	ADS	16 mm	Green
MELAG	63500	16 mm	Green
GETINGE	502887400	16 mm	Black

- ③ NSK handpieces with this symbol are compatible with medical washer-disinfectors.  VA2-LUX-HP Ti, VA2-LUX-HP, VA2-HP can be cleaned and disinfected with a washer-disinfector.

Due to the variation in cleaning/disinfection cycles and cleaning agents available from the various manufacturers, NSK has validated the Getinge 46 Series Washer Disinfector Cycle P7-Anesthesia for automatic processing/reprocessing of the product.

Only validated procedures must be used for cleaning and sterilization. When using reprocessing procedures that are different from those described in this manual, those procedures must be validated by the respective practice or hospital using a washer-disinfector that conforms to ISO 15883-1.

NOTICE

- All water used in cleaning is tap water of drinking water quality unless otherwise indicated.

Perform the main wash under the following conditions. Use a Getinge 46 Series washer-disinfector (or equivalent device):

Pre-Wash 1: Cold water (<43 °C (110 °F)) for 2 min

Pre-Wash 2: None - Not Applicable

Enzyme Wash: Hot tap water and 60 ml of heated (35 °C (95 °F)) neutral medical enzymatic detergent for 3 min at 49 °C (120 °F)

Detergent Wash: None - Not Applicable

Rinse 1: Hot tap water for 2 min

Rinse 2: Hot tap water for 2 min

RO Final Rinse: Hot tap water for 32 min at 75 °C (167 °F)

Drying: 20 min

Upon removal from the washer-disinfector, if residual moisture is present, wipe off with dry cloth or blow it off with compressed air (≤ 0.35 MPa) until there is no moisture in the interior and exterior.

- 4 Under appropriate lighting (500 lx or more), check that there is no dirt. If there is dirt, repeat from step 1 until the dirt is removed.

Proceed to sterilization. $\Rightarrow$ 5-6 Packaging, Sterilizing, and Drying

Figure 30

⚠ CAUTION

- Handpieces must be removed from the washer-disinfector immediately (within 1 hour) after the cleaning and drying cycle is complete to prevent corrosion.
- After washing with washer-disinfector, dry the product until all internal moisture is completely removed.
- If any washer-disinfector moisture remains inside the product, it could reduce the effect of lubrication and could cause corrosion inside the product.

5-5 Cleaning glass rods

If dust or debris adheres to the glass rod, carefully wipe it off using a cotton swab dampened with disinfectant alcohol.


Figure 31

⚠ CAUTION

- When cleaning the end face of the glass rod, do not use needles or sharp tools. Scratches can reduce light transmittance. If it becomes dark due to scratches, please contact an authorized NSK dealer.

Type	Gravity Displacement	Pre-Vacuum (Dynamic Air Removal)
Temperature	132 °C (270 °F)	132 °C (270 °F)
Holding Time	15 min	4 min
Drying Time	30 min	30 min

5-6 Packaging, Sterilizing, and Drying

- 1 Insert the handpiece into an FDA-approved sterilization pouch (es) that conforms to ISO 11607-1 and seal the pouch (es). Perform steam sterilization with the following conditions

CAUTION

- Use an FDA-approved steam sterilizer to perform sterilization.
- Follow local rules, regulations, and guidelines regarding the reprocessing of devices.
- Remove the product from the sterilizer immediately after sterilization is complete. Be careful not to burn yourself, because the handpiece will be very hot immediately after sterilization.
- Place each product in an individual sterilization pouch for steam sterilization. This is to prevent possible discoloration and damage to the product due to chemical residue from other instruments.

CAUTION

- Do not use sterilization pouches containing water-soluble adhesive components such as PVA (polyvinyl alcohol). Adhesive components washed out during sterilization may enter the product, causing problems such as poor rotation, sticking and failure to operate correctly.
- Clean the handpiece prior to sterilization. If blood remains on the internal surface, it can become clotted and cause product failure.
- Do not heat or cool the product too quickly. A rapid change in temperature could cause damage to the product.

CAUTION

- Be sure to use sterilizers that can perform sterilization up to 135 °C (275 °F). In some sterilizers, the chamber temperature may exceed 135 °C (275 °F). Do not use these sterilizers because failure of the handpiece could occur. Contact the sterilizer manufacturer for detailed information about cycle temperatures.
- Steam sterilization is recommended for the product. The validity of other sterilization methods (such as plasma sterilization or EOG sterilization) is not confirmed.

NOTICE

- Repeated steam sterilization may cause the handpiece to discolor due to heat. This is a characteristic of the product and does not affect its quality.

5-7 Storage

Store the product in a dry, clean location.

CAUTION

- Store the product in a well ventilated place out of direct sunlight and within the range of temperature, humidity and pressure specified in “8 Specification”.
- Sterilization is not guaranteed after the sterilization retention period specified by the manufacturer and seller of the sterilization pouch has elapsed. If the sterilization retention period has elapsed, perform sterilization again with a new sterilization pouch.

CAUTION

- Store the product in a place where it will not be adversely affected by air pressure, temperature, humidity, ventilation, sunlight, dust, salt, or sulfur-containing air.

6 Periodical Maintenance Checks

Perform periodical maintenance checks every three months, while referring to the check sheet below. If any abnormalities are found, contact an authorized NSK dealer.

Product Appearance	Check for any scratches, damage, or abnormalities.
Connection parts	Check each connection for scratches, damage, and deformation.
Connection condition	Connect the handpiece cord and tip, and check that there is no looseness in the connections.
Vibration	After connecting each part, press the foot control and let the handpiece and tip vibrate for about 1 minute, then check the following items: – Water is being supplied. – The light turns on (for handpieces with lights). – There are no abnormal vibrations, noises, or overheating. – The vibration stops normally after releasing the foot control.

7 Troubleshooting

*Cannot be fixed at dental clinics. Be sure to contact an authorized NSK dealer.

Problem	Cause	Remedy
There is no vibration or the vibration is poor.	The threads of the ultrasonic handpiece and tip are not clean.	Wipe them clean.
	The threads of the ultrasonic handpiece are worn out as a result of long term use.	Contact an authorized NSK dealer*.
	Failure of the vibration device in the ultrasonic handpiece	

Problem	Cause	Remedy
The ultrasonic handpiece produces a loud noise.	Power has not been correctly adjusted for the tip.	Adjust the power so that it does not exceed the set upper limit value mentioned in the supplied operation manual of the tip. ⇒ 1-2 Related manuals
	The tip is not tightened firmly.	Re-tighten the tip with the tip wrench by rotating the tip wrench until it makes a click sound and cannot rotate any further. If the problem still persists, replace the tip wrench as it may be worn out.
	Failure of the ultrasonic handpiece's vibration device or control unit's internal components.	Contact an authorized NSK dealer*.

Problem	Cause	Remedy
The ultrasonic handpiece is overheating.	Power has not been correctly adjusted for the tip.	Adjust the power so that it does not exceed the set upper limit value mentioned in the supplied operation manual of the tip. ⇒ 1-2 Related manuals
	The tip is not tightened firmly.	Re-tighten the tip with the tip wrench by rotating the tip wrench until it makes a click sound and cannot rotate any further. If the problem still persists, replace the tip wrench as it may be worn out.
	The threads of the ultrasonic handpiece and tip are not clean.	Wipe them clean.
	The threads of the ultrasonic handpiece are worn out as a result of long term use.	Contact an authorized NSK dealer*.
	Failure of the ultrasonic handpiece's vibration device or control unit's internal components.	

Problem	Cause	Remedy
Damage to the appearance or metal parts of the ultrasonic handpiece.	The ultrasonic handpiece has been dropped or deteriorated, or unspecified chemicals have been used.	Contact an authorized NSK dealer*.

8 Specification

8-1 Handpiece specification

Model	VA2-LUX-HP Ti	VA2-LUX-HP	VA2-HP
Lighting	Glass rod	Glass rod	-
Vibration frequency	28 - 32 kHz	28 - 32 kHz	28 - 32 kHz
Dimensions	Ø20 mm×103 mm	Ø20 mm×103 mm	Ø17.6 mm×100 mm
Weight	58.5 g	53.5 g	42.5 g

	Temperature	Humidity	Atmospheric pressure
Use Environment	+ 10 - 35 °C (50-95 °F)	30 - 75 %	700 - 1060 hPa
Store Environment	- 10 - 50 °C (14-122 °F)	10 - 85 %	500 - 1060 hPa

Use outside of these limits may cause malfunction.

56 Avoid any risk of freezing of irrigation solution. When operated outside the above range, there is a risk of failure.

8-2 Classification of equipment

Type of protection against electric shock: Refer to the operation manual of the ultrasonic system.

- Degree of protection against electric shock: Refer to the operation manual of the ultrasonic system.
- Method of sterilization or disinfection recommended by the manufacturer:
 - Autoclave sterilization
- Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide:
 - Equipment NOT suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide
- Mode of operation: Refer to the operation manual of the ultrasonic system.

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NSK America Corp. www.nskdental.com